

BARIATRIC BEDPAN

H133-15

The Bariatric Bedpan by Medegen Medical Products is an enhanced fracture bedpan design with reinforced construction and embedded handles that provides a comfortable surface for the patient. The unique color makes it easy for the clinician to readily identify the bariatric bedpan.



- Low profile to improve usability
- Integrated handles for nurse safety and reduced injuries
- Center support for improved patient weight bearing
- Broad flat rim for improved patient support and comfort
- Holds 2 qt. of waste
- Easy to clean and sanitize with integrated pour spout
- Dimensions: 15.75" L x 14.75" W x 2.75" H (overall); 0.5" H (back rim)
- Weight: 13.3 oz.
- Weight Capacity: 750 lbs. dead weight (Patient weight approximately 1200 lbs.)
- Can be used for patients as low as 200 lbs.
- Recyclable 

To view our products online, please visit us at www.medegenmed.com or directly with your device by scanning the QR code to the right.



**MEDEGEN**
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

Customer Service Fax 631-714-6311

E-mail cservice@medegenmed.com