

INTRODUCING

the **KALIKART™** series of
HEALTHCARE CARTS

Quiet, Robust, Easy to Move | Designed to Meet All of Your Needs



Carts are used in just about every area of Healthcare whether it be Hospitals, Non-Acute Care facilities, Surgery Centers, Doctor's offices or in a Lab.

Carts are essential to Healthcare providers allowing the transportation and storing of medical supplies safely and conveniently.

Medegen Medical KaliKart offerings are multi-functional and can be configured to meet your specific needs. KaliKart products are constructed with maintenance free ABS plastic, aluminum and powder coated steel.

Which means/results in:

- Rust Proof
- Longer life and lower overall cost
- Quieter opening and closing providing a better Patient Environment
- Easier to clean and disinfect
- Easy to maintain

KaliKart units also come with:

- 5" casters – locking and non-locking
- Integrated dividable draws
- Locking system

KaliKart work surfaces come in various heights:

33", 37", 40" and 44".

We offer anti-microbial color coded draw pulls.

Color options are:

- | | | |
|---------|----------|----------|
| ● Red | ● Blue | ● Yellow |
| ● Green | ● Orange | ● Purple |
| ● Beige | ● Grey | ○ White |

We have a full line of accessories!





KALIKART™ ANESTHESIA SOLUTION

- Top Surface L23.6" x W17.7"
- Slide out Working Shelf
- 4 Casters: 2 locking/2 non-locking
- 2 3.2"H Drawers with DIY Dividable Set
- 2 6.3"H Drawers with DIY Dividable Set
- 1 9.5"H Drawer with DIY Dividable Set
- Monodose Cassette (S)
- Bridge UD Dose shelf
- 8708 5.4 Sharps container w/bracket
- Side Cabinet with 4 Drawers
- DB Triple Waste bin
- Key Lock for drawer

KALIKART™ EMERGENCY SOLUTION

- | | |
|---|---|
| <ul style="list-style-type: none"> ■ Top Surface L23.6" x W17.7" ■ Slide out Working Shelf ■ 4 Casters: 2 locking/2 non-locking ■ 2 3.2"H Drawers with DIY Dividable Set ■ 2 6.3"H Drawers with DIY Dividable Set ■ 1 9.5"H Drawer with DIY Dividable Set | <ul style="list-style-type: none"> ■ IV Pole ■ Defibrillator Shelf ■ Oxygen Tank Holder ■ Seal Lock |
|---|---|



KALIKART™ MEDICATION SOLUTION

- Top Surface L23.6" x W17.7"
- Slide out Working Shelf
- 4 Casters: 2 locking/2 non-locking
- 1 9.5"H Drawer with DIY Dividable Set
- 6 Dose Set
- Front Push Handle
- Key Lock

KALIKART™ TREATMENT SOLUTION

- Top Surface L23.6" x W17.7"
- Slide out Working Shelf
- 4 Casters: 2 locking/2 non-locking
- 2 3.2"H Drawers with DIY Dividable Set
- 2 6.3"H Drawers with DIY Dividable Set
- 1 9.5"H Drawer with DIY Dividable Set
- Key Lock



KALIKART™ WORKSTATION SOLUTION

- Top Surface L23.6" x W17.7"
- Slide out Working Shelf
- 4 Casters: 2 locking/2 non-locking
- 2 3.2"H Drawers with DIY Dividable Set
- 2 6.3"H Drawers with DIY Dividable Set
- 1 9.5"H Drawer with DIY Dividable Set
- PC Arm and Holder
- Front Push Handle
- Key Lock



KALIKART™ MED-SURG SOLUTION

- Top Surface L23.6" x W17.7"
- Slide out Working Shelf
- 4 Casters: 2 locking/2 non-locking
- 2 3.2"H Drawers with DIY Dividable Set
- 2 6.3"H Drawers with DIY Dividable Set
- 1 9.5"H Drawer with DIY Dividable Set
- Side Cabinet with 4 Drawers
- Double Waste Bin
- IV Pole
- Push Ball Handle
- Key Lock



KALIKART™ ACCESSORIES



O₂ Holder/Side Cabinet Set



Side Working Tray



Side Tray with Compartments



Defibrillator Shelf




Dividable Drawers



Trash & Biohazard Waste Receptacles



CPR Massage Board




Corner Push Handle




RFID Sensor Lock with Corner Push Handle



Short IV Pole




Basket



Key Lock



Storage Buckets



Seal Lock



Locking Casters



Monodose Cassettes

Call 800.511.6298 or on the web: www.medegenmed.com

OTHER GREAT PRODUCTS FROM MEDEGEN


HOSPITAL FURNITURE



*See our full line of furniture in this
informative product brochure!*



CHEVALIER Hydraulic STRETCHER



*See all the features in
our comprehensive trifold
product brochure!*



*To view our products online, please visit us at
www.medegenmed.com or directly with your
device by scanning the QR code to the right.*



KaliKart – Medegen Medical Products offer 1 year warranty for ABS plastics cart parts and for our accessories/components and a 3 year warranty for the carts main frame/metal parts.



360 Motor Pkwy, Suite 800 • Hauppauge, NY 11788 • 800.511.6298
Customer Service Fax: 631.714.6311 • Customer Service: CSservice@medegenmed.com
International Customer Service: IntCSservice@medegenmed.com

www.medegenmed.com

© 2022 Medegen Medical Products Lit# 0080D

