



HOSPITAL FURNITURE

*Quality, Comfort and
Convenience*



HOSPITAL FURNITURE

RECLINER-SLEEPER

1 Operating Handles

Operating Handles are on the left and right side of the recliner, making it easy to move or assisting with converting the recliner to a sleeper.



2 Safety Knobs

Safety Knobs are located on the left and right sides of the recliner-sleeper. Either knob can easily convert the recliner to a sleeper and back to a recliner.




3 Sleeper

Recliner-Sleeper has an ergonomic sleeping surface with drop down armrests providing an unobstructed sleeping area.



4 Gas Spring Back

The recliner has a Gas Spring back making the transition from recliner to sleeping surface effortless.



5 Arm Rests

Arm rests are made of soft comfortable foam.



6 Braking Caster

Recliner comes with 4 x 3" casters, 2 are locking and are toe operated.



Medegen Recliner-Sleepers are offered in 3 different sizes and 2 colors:

SMALL

Chair: L 34" x W 24" x H 39"
Bed: L 77" x W 24" x H 18"
Sleeping surface: L 73" x W 21"

MEDIUM

Chair: L 34" x W 27" x H 39"
Bed: L 79" x W 27" x H 18"
Sleeping surface: L 75" x W 22"

LARGE

Chair: L 35" x W 30" x H 42"
Bed: L 85" x W 30" x H 18"
Sleeping surface: L 80" x W 24"

COFFEE

BEIGE

PATIENT BEDSIDE CABINETS

Our patient Bedside Cabinets come assembled and are offered in wood grain laminate or ABS plastic.

- Color offered in wood laminate: Maple
- Color offered in ABS: Blue
- ABS plastic Bedside Cabinet top has a lip on 3 sides preventing item from falling off
- Left or Right mountable hanger for either Bedside Cabinet
- Wood laminate comes with a slide out shelf, draw, closet and a vented drawer
- ABS plastic comes with a slide out shelf, drawer, closet and a vented shelf
- 4 x 2" casters, front casters are locking
- Dimension: L 22" x W 19" x H 30"



PATIENT OVERBED TABLE

Medegen's Overbed Tables are designed to be:

- STABLE – with an “H” shaped base and four casters (2 x 1" and 2 x 2")
- DURABLE – we offer two options: wood grain laminate or ABS Plastic
- FUNCTIONAL – our Overbed Table has a gas spring lifting column with a fingertip self-locking handle enabling Healthcare workers to adjust the table height from 30" up to 41". Accommodating the Patient Bed, Stretcher, Chair/Recliner or a Wheelchair
- Overbed Table Dimensions are: L 35" x W 17"
- Overbed Tables are available in ABS plastic color Blue and Wood Laminate color Light Maple
- Table Weight capacity is 60 LBS



OTHER GREAT PRODUCTS FROM MEDEGEN


the **KALIKART™** series of **HEALTHCARE CARTS**



*See our full line of carts in this
informative trifold product brochure!*



CHEVALIER **Hydraulic STRETCHER**



*See all the features in our
comprehensive trifold
product brochure!*



*To view our products online, please visit us at
www.medegenmed.com or directly with your
device by scanning the QR code to the right.*

