

STACKABLE SHARPS CONTAINERS

5.4 Qt. **SHARPS-tainer®**

FEATURES/BENEFITS:

- Puncture and Leak Resistant
- Locking Lid
- Autoclave/Incinerator Safe up to a max of 280°F (137°C)
- Fits Industry 5 and 5.4 qt. cabinets
- Increases 5 qt. cabinet capacity to 5.4 qt.
- *Typical areas of use:*
OR, ER, Surg. Centers, Patient Rooms, Phlebotomy Labs, Chemistry, Hematology, Pathology, Toxicology, Microbiology, Immunization, Drug Screening Labs and wherever Biohazardous Sharps exists.

NEW!



8708



8708C



8708TY



8708T



To view our products online, please visit us at www.medegenmed.com or directly with your device by scanning the QR code to the left.

MEDEGEN
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036
Toll Free 800-511-6298 • Customer Service Fax 631-714-6311
E-mail cservice@medegenmed.com



Stackable Features & Benefits **SHARPS-tainer®**

- Stackable design provides a stable platform during use, economical storage and lower shipping costs.
- Wide horizontal drop lid models accommodates larger syringes
- Tortuous path lids help limit hand exposure for added security
- Transparent lids facilitates quick visual identification of fill level
- Translucent containers facilitate easy side viewing of content level
- Compact containers are convent for easy transport
- Large containers accommodate endoscopic accessories, elongated disposable instruments and probes
- Chemo Waste Containers (WHITE) are CYTA labeled for easy identification

Code	Container/Print Color	Size (inches)	Capacity	Qty/cs
A • 8716	Red/Black	13.5 x 13.5 x 25.5	16gal.	6
B • 8708TY	Yellow/Black	10.5 x 3.5 x 10.5	5.4qt.	12
C • 8708	Red/Black	10.5 x 3.5 x 10.5	5.4qt.	12
D • 9752	White/Black & Yellow	13.5 x 13.5 x 25.5	16gal.	6
E • 8710	Red/Black	12.25 x 6 x 13.5	3gal.	12
F • 8704T	Translucent Red/Black	10 x 7 x 9.5	2gal.	24
G • 8705TY	Yellow/Black	13.75 x 13.75 x 14	8gal.	10
H • 8710TF	Red/Black	12.5 x 6 x 12	3gal.	12
I • 8704TY	Yellow/Black	10 x 7 x 9.5	2gal.	24
J • 8702T	Translucent Red/Black	3.5 x 3.5 x 7	1qt.	72
K • 8702	Red/Black	3.5 x 3.5 x 7	1qt.	72
L • 8702TY	Yellow/Black	3.5 x 3.5 x 7	1qt.	72
M • 8705T	Translucent Red/Black	13.75 x 13.75 x 14	8gal.	10

Code	Container/Print Color	Size (inches)	Capacity	Qty/cs
N • 8707	Red/Black	10 x 7 x 11.25	2gal.	10
O • 8706	Red/Black	10 x 7 x 6.75	1gal.	12
P • 8705	Red/Black	13.75 x 13.75 x 14	8gal.	10
Q • 8707TY	Yellow/Black	10 x 7 x 11.25	2gal.	10
R • 8706T	Translucent Red/Black	10 x 7 x 6.75	1gal.	12
S • 9751	White/Black & Yellow	13.75 x 13.75 x 14	8gal.	10
T • 8704	Red/Black	10 x 7 x 9.5	2gal.	24
U • 8707T	Translucent Red/Black	10 x 7 x 11.25	2gal.	10
V • 8703T	Translucent Red/Black	10 x 7 x 5	1gal.	24
W • 8708T	Red/Black	10.5 x 3.5 x 10.5	5.4qt.	12
X • 8703	Red/Black	10 x 7 x 5	1gal.	24
Y • 8708C	Clear/Black	10.5 x 3.5 x 10.5	5.4qt.	12



Non-Stackable Features & Benefits **SHARPS-tainer®**

- Funneled vertical drop and extended neck design helps minimize needle exposure and over filling
- Stable base during use
- Extended handles keep fingers away from opening during transport
- Needle removal port on 158 and 178 which provides convenient and compact disposal of sharps
- Translucent sharps have wrap around window view label which provides easy visualization of fill level
- Chemo Waste models (white) are labeled with CYTA for easy identification

Code	Description	Container/Print Color	Size (inches)	Capacity	Qty/cs
A • 195	Tall tray Size, Locking Screw Cap	Clear/Orange & Black	5 x 4 x 10	2qt.	20
B • 182	Medium Size, Lacking Cap	Red/Orange & Black	10 x 6 x 11	6.2qt.	12
C • 180S	Large Size, Locking Screw Cap	Clear/Orange & Black	10 x 6 x 15	10qt.	12
D • 932	Large Size, Large Opening, Locking Cap	White/Orange & Black	12 x 9 x 18	23qt.	4
E • 931W	Standard Size, Large Opening, Locking Cap	White/Orange & Black	12 x 9 x 11	12qt.	12
F • 158	Short Tray Size w/Needle Remover	Red/Orange & Black	3.5 x 3.5 x 5	.7qt.	40
G • 185S	Short Tray Size, Locking Cap	Red/Orange & Black	3.5 x 3.5 x 5	.7qt.	40
H • 199	Short Tray Size, Locking Screw Cap	Clear/Orange & Black	5 x 4 x 7.5	.7qt.	40
I • 184	Small Size, Locking Cap	Red/Orange & Black	10 x 6 x 9	4.8qt.	12
J • 187	Extra Large Wide Opening, Locking Cap	Red/Orange & Black	12 x 9 x 18	23qt.	4
K • 186S	Large Wide Opening, Locking Cap	Clear/Orange & Black	12 x 9 x 11	12qt.	12
L • 184R	Rectangular, Locking Cap	Red/Orange & Black	10.25 x 6 x 7.75	4.7qt.	12
M • 187S	Extra Large Wide Opening, Locking Cap	Clear/Orange & Black	12 x 9 x 18	23qt.	4
N • 186W	Large Wide Opening, Locking Cap	Red/Orange & Black	12 x 9 x 11	12qt.	12
O • 184S	Small Size, Locking Screw Cap	Clear/Orange & Black	10 x 6 x 9	4.8qt.	12
P • 180	Large Size, Locking Cap	Red/Orange & Black	10 x 6 x 15	10qt.	12
Q • 178	Tall Tray Size, Locking Cap	Red/Orange & Black	3.5 x 3.5 x 10	1.7qt.	20
R • 930	Small Size, Locking Cap	White/Orange & Black	3.5 x 3.5 x 10	1.7qt.	12
S • 185	Tall Tray Size w/Needle Remover	Red/Orange & Black	3.5 x 3.5 x 10	1.7qt.	20
T • 182S	Medium Size, Locking Screw Cap	Clear/Orange & Black	10 x 6 x 11	6.2qt.	12



INTEPLAST GROUP
Healthcare