

TRI-HOLD BEDPAN

H113 Series

*A stackable, convenient
bedpan solution*

FEATURES/BENEFITS:

- Textured seat for patient comfort and convenience
- Designed to fit over the commode/toilet for collection
- Pour areas for easy emptying
- Three integrated handles for easy placement/removal
- Maximizes storage options, units nest together, stackable
- Packaging:
50/case, bulk
- Available in four colors: Gold (-05),
Turquoise (-07),
Dusty Rose (-10),
and Gray (-11)
- Recyclable 
- Made in the USA




H113-05

H113-07

H113-10

H113-11

 **MEDEGEN**
Medical Products
An Inteplast Group Company


209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298 • Customer Service Fax 631-714-6311

E-mail cservice@medegenmed.com



To view our products online, please visit us at
www.medegenmed.com or directly with your
device by scanning the QR code to the left.



A H102
B H100
C 00080
D H101
E H114
F 89020
G H111
H H103
I 00083



A 00081
B 00085
C H120
D H112
E H133
F H7870
G H113
H 89010
I H131
J H122

Bedpans (Single Patient Use)

Code#	Description	Capacity	Weight Capacity‡	Color	Qty/cs	Packaging
H100	Fractured	34 oz.	250lbs	Various	50	Bulk
H100BG-05	Fractured	34 oz.	250lbs	Gold	50	Individually Bagged
H101	Lg. Fractured	2.5L	350lbs	Various	12	Bulk
H102	Fractured	34 oz.	250lbs	Various	12	Bulk
H103-10	Fracture, Deluxe	2 qt./1.9L	350lbs	Dusty Rose	50	Bulk
H111	Stackable w/Skirt	1.75 qt./1.66L	250lbs	Various	50	Bulk
H112	Stackable	2 qt./1.9L	250lbs	Various	50	Bulk
H113	Commode, Stackable w/Handles	2 qt./1.9L	300lbs	Various	50	Bulk
H114-05	Stackable w/Lg. Skirt	2.25 qt./2.13L	300lbs	Gold	25	Individually Bagged
H120	Pontoon	2 qt./1.9L	250lbs	Various	20	Bulk
H120BG	Pontoon	2 qt./1.9L	250lbs	Various	20	Individually Bagged
H122	Pontoon w/Cover	2 qt./1.9L	250lbs	Various	10	Individually Bagged
H130-05	Conventional	3 qt./2.84L	350lbs	Gold	12	Individually Bagged
H131	Conventional	3 qt./2.84L	350lbs	Various	12	Bulk
H131BG-10	Conventional	3 qt./2.84L	350lbs	Dusty Rose	12	Individually Bagged
H133-15	Bariatric	2 qt.	750lbs	Orange	6	Bulk
H7870	Commode, Stackable w/Handles	2 qt./1.9L	300lbs	Various	50	Bulk

Bedpans (Sterilizable Options)

Code#	Description	Capacity	Weight Capacity‡	Color	Qty/cs	Packaging
00080	Conventional	2.25 qt./2.25L	300-350lbs	Blue	6	Bulk
00081	Fracture	1 qt.	250lbs	Blue	6	Bulk
00083	Lg Fracture	2.5L	350lbs	Blue	6	Bulk
00085	Pontoon/Contour	2 qt./1.9L	250lbs	Blue	12	Bulk
89010	Bedpan	3 qt./2.84L	500lbs	Stainless Steel	1	Individual
89020	Fractured	34 oz.	500lbs	Stainless Steel	1	Individual

*Sterilizable products can withstand temperature up to 275°F (135°C), steam sterilization ‡ Weight Capacity is "deadweight" and not patient weight. Bedpans hold the weight of the patient buttocks and mid-trunk of their body. The head, arms, legs and chest are supported by the bedding.

FEATURES/ BENEFITS:

- Single Patient Use and/or Reusable
- Multiple color options
- Sterilizable options
- Packaging Options
 - Bulk
 - Individually bagged
 - Stainless Steel - Individual



INTEPLAST GROUP
Healthcare

© 2016 Medegen Medical Products
LIT# 0047B