

NEW!

DUAL TEST AUTOCLAVE BAGS

New Impact (Dart) and Tear (Elmendorf) compliant Autoclave Bags provide trusted high heat performance of traditional Medegen Autoclave bags with the added peace of mind that comes with Dual-Tested compliance.

Autoclaving is a very common way to sterilize infectious waste through the use of high temperature, high pressure steam over a specific time period. Sterilized waste is less expensive to dispose of than infectious biohazard waste. New Medegen Dual-Tested Autoclave Bags are ASTM D1922 Elmendorf Tear tested to exceed 480g of tear strength (MD and TD) and ASTM D1709 tested to exceed 165g dart drop strength. These bags comply with US DOT 49CFR Sections 173-197 and specific state regulations.

FEATURES/BENEFITS:

- Dual Tested compliant
- Operating temperature tolerance up to 285°F
- Integral steam processing indicator that turns brown when autoclaved
- Available in various sizes and 2 film colors (Red or Clear)
- Customization available – Contact your local Medegen Medical Products Sales Person.



To view our products online, please visit us at www.medegenmed.com or directly with your device by scanning the QR code to the left.



Item No.	Size (in)	Gauge	Film Color	Print	Bags/Cs
AC1419C	14x19	2 mil	Clear	Black	500
AC1419R	14x19	2 mil	Red	Black	500
AC1923C	19x23	2 mil	Clear	Black	200
AC1923R	19x23	2 mil	Red	Black	200
AC2530C	25x30	2 mil	Clear	Black	200
AC2530R	25x30	2 mil	Red	Black	200
AC2535C	25x35	2 mil	Clear	Black	200
AC2535R	25x35	2 mil	Red	Black	200
AC3138C	31x38	2 mil	Clear	Black	100
AC3138R	31x38	2 mil	Red	Black	100
AC3847C	38x47	2 mil	Clear	Black	100
AC3847R	38x47	2 mil	Red	Black	100

MEDEGEN
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036
Toll Free 800-511-6298 • Customer Service Fax 631-714-6311
E-mail cservice@medegenmed.com