

# FEMALE URINAL


**H145-01**

*Supports your initiative to reduce  
Urinary Tract Infections (UTI)*

## A Better Alternative to a Fracture Bedpan!

- More comfortable for the patient
- Patient can use at their own convenience
- Nurse time saver
- Less prone to spillage
- Comfortable, elongated body design

## FEATURES/BENEFITS:

- Contoured smooth opening
- Contents easily measured
- Specially designed handle conveniently hangs on most beds or with H143-01 Uri-Clip™
- Pronounced neck angle helps confine contents when urinal is positioned horizontally
- Graduated: 25 cc and 1 oz. increments up to 1000 cc/32 oz
- Packaging: 36/case, bulk
- Recyclable 



 **MEDEGEN**  
Medical Products  
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298 • Customer Service Fax 631-714-6311

E-mail [cservice@medegenmed.com](mailto:cservice@medegenmed.com)



To view our products online, please visit us at  
[www.medegenmed.com](http://www.medegenmed.com) or directly with your  
device by scanning the QR code to the left.

# Your Source For URINALS



## FEATURES/ BENEFITS:

- Single Patient Use and/or Reusable
- Sterilizable options
- Translucent urinals for content visualization and measurement ease
- Wide base design for stability empty or full
- Odor shield (lid) options

### Packaging Options:

- Bulk
- Individually bagged

### Medegen Medical Products

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

E-mail [cservice@medegenmed.com](mailto:cservice@medegenmed.com)

Web [www.medegenmed.com](http://www.medegenmed.com)

### Urinals - Male (Single Patient Use)

Code#	Color	Capacity	Qty/cs	Packaging
H140-01	Translucent	1qt/1000cc	50	Bulk
H140CR-01	Translucent w/Comfort Ring	1qt/1000cc	50	Bulk
H141-01	Translucent	1qt/1000cc	50	Individually Bagged

### Urinals - Male (Sterilizable\*)

Code#	Color	Capacity	Qty/cs	Packaging
00090	Blue	1.25qt/1180cc	12	Bulk w/Liner
00095	Blue	1qt/1000cc	12	Individually Bagged
00195	Translucent w/Blue Cover	1qt/1000cc	12	Individually Bagged
89150	Stainless Steel	7/8qt / .828L	1	Sold Individually

\*Sterilizable products can withstand temperature up to 275°F (135°C), steam sterilization

### Urinals - Female (Single Patient Use)

Code#	Color	Capacity	Qty/cs	Packaging
H145-01	Translucent	1qt/1000cc	36	Bulk

### Uri-Clip™ Urinal Holder

Code#	Color	Capacity	Qty/cs	Packaging
H143-01	White	N/A	48	Bulk