

In the medical profession, every procedure is critical. We dedicate every working day to the design, manufacturing, packaging and distribution of medical and surgical disposable products using stringent quality controls. Through this meticulous attention to detail, we deliver healthcare products with the assurance that each one will be crafted to an exacting standard of quality. Every item we offer is designed for superior performance, patient comfort and safety.

Patient Bedside Products

- Bedpans
- Central Supply - Stainless Steel
- Cleansing & Irrigating Sets
- Denture Cups
- Emesis Basins
- Emesis Bags
- Irrigating & Suctioning Syringes
- Medicine Cups
- Mobility Aids
- Pitchers and Carafes
- Procedure Trays
- Specimen Collectors
- Stainless Steel Essentials
- Triangular Graduates
- Urinals
- Wash Basins

Specimen Collection and Laboratory Products

- 24-Hour Collection Containers
- ClikSeal™ Specimen Containers
- Graduated Measures
- Histology / Pathology Containers
- Intake and Measuring Containers
- Laboratory Containers & Bottles
- Midstream Catch Kits
- Non-Sterile Specimen Containers
- Other Laboratory Products
- Patient Specimen Collectors
- Snap Cap Specimen Containers
- Spouted Collection Containers
- Sterile Specimen Containers
- Transport Swabs
- Tri-Pour Beakers & Lids
- Urinary Drainage Bags

Collection and Containment

- Biohazardous Waste Containment
- Autoclavable Biohazardous Waste Bags
- Biohazardous Waste Bags
- Non-Stackable Sharps Containers
- Stackable Sharps Containers
- Fluid and Spill Management
- Clean-Up Kits
- Solidifiers
- Institutional Trash Containment
- High Density Trash Can Liners
- Low Density Trash Can Liners
- Laundry and Linen Containment
- Hardware and Accessories
- Laundry and Linen Collection Bags
- Water Soluble Laundry Bags
- Specialty Products
- Lab Specimen Transport Bags
- Medi-Waste® System
- Medical Equipment and Cart Dust Covers
- Other Collection and Containment
- Essentials
- Patient Belongings Bags
- Sterility Maintenance Covers
- Zip Closure Bags




For detailed information about any of our products or services, contact your Medegen customer service representative at **1.800.511.6298**. To view our products online, please visit us at **www.medegenmed.com** or directly with your device by scanning the QR code to the left.



Medegen Medical Products
209 Medegen Drive
Gallaway, TN 38036

CUSTOMER SERVICE 800.511.6298
EMAIL CService@medegenmed.com
FAX 631.714.6311
INTERNATIONAL IntCService@medegenmed.com
WEB www.medegenmed.com



Lightweight | Adjustable | Comfortable | Secure

Medegen Medical Product’s mobility aids are designed to help continue patient’s active lifestyle and improve mobility. Our crutches, canes and walkers assist in maintaining upright ambulation by providing: improved stability, reduced lower-limb loading and assist with generating movement.

walkers

- Deep, wide frame allows for patient’s full stride and mobility
- 4 bottom rubber tips adds a safe grip to prevent falling or sliding
- Front cross bar and lower side brace add extra support and stability
- Easy push button adjusts height smoothly in 1" increments, offering a personalized fit that accommodates a broad range of patients’ heights
- Compact design folds to a 4" depth for easy storage and transport



Part #	Description	Weight Capacity
WLK1001	Walker: Folding 1-Button 32-36	300 lbs
WLK2001	Walker: Folding 2-Button 32-39	500 lbs

Case Qty: 4. Packaged individually in transparent bag. Sizing instructions included. No assembly required.

Use your mobile device to view how to fit crutches, canes and walkers properly by scanning the QR code to the right.

www.medegenmed.com/docs/mobility-sizing.pdf



WLK1001 Height: 32"-36"

- 1-button audible click folding with covered button to protect skin from piercing
- Easy one hand push button for quick and simple folding & unfolding
- Audible click provides an added safety feature to ensure walker is either completely opened or completely closed



DELUXE
500 lb weight capacity

WLK2001 Height: 32"-39"

- 2-button audible click folding mechanisms allow patient to close one side at a time for maneuvering through tight spaces and corners
- The locking audible click mechanism locks one side to ensure sturdiness and safety as weight is distributed to one side
- 500 lbs weight capacity supports larger patients

crutches

300 lb weight capacity

- Anodized aluminum construction: lightweight, durable and easy to clean
- Wide range of sizes to accommodate all patients’ heights; 4 sizes from Pediatric to Tall
- Adjustable push button with convenient printed 1" increment height scale, ensures a quick & easy customization of height and accurate fitting
- Rubber armrest pad and hand grip offer maximum comfort and the best possible pressure distribution

Part #	Description	Weight Capacity	Patient Height
CRUP001	Crutch: Aluminum Pediatric 31.5-37.5	300 lbs	4'0" - 4'6"
CRUY002	Crutch: Aluminum Youth 37-45	300 lbs	4'6" - 5'2"
CRUA003	Crutch: Aluminum Adult 45-53	300 lbs	5'2" - 5'10"
CRUT004	Crutch: Aluminum Adult Tall 53-61	300 lbs	5'10" - 6'6"



Case Qty: 8. Shrink wrapped in pairs. Sizing instructions included. No assembly required.

canes

- Height adjusts in 1" increments from 30" to 39" for a convenient personalized fit
- Simple push button height with locking collet
- Slip resistant bottom rubber tip protects from sliding and slipping
- Contoured PVC handgrip offers a comfortable and secure hand grip
- Easy and quick disassemble is perfect for storing and traveling

Part #	Description	Weight Capacity
CNE0001	Cane: Aluminum Adjustable 30-39	300 lbs

Case Qty: 6. Packaged individually in transparent bag. Sizing instructions included. No assembly required.

