

NEW!

Super Absorbent PADS

for **ANY** Specimen Bag

Minigrip 3"x5" Super Absorbent Pads have always been available pre-stuffed into Minigrip specimen bag items SBL2AP69B and UF95-BDES. These super-absorbent pads are now available as an item in both the Minigrip and Medegen product lines.

BAG FLEXIBILITY

- Absorbent pads can be placed into ANY specimen bag.
- More than one pad can be inserted into a specimen bag for extra absorbency.

CONVENIENCE PACKAGING

- Absorbent pads are packed 100 / pack, 10 packs / case – just like many Medegen, Minigrip and competitor specimen bags. User can use one pack of pads with one pack of bags.

SUPER ABSORBENCY

- This compact 15 square inch pad delivers 25 – 35ml of urine and 100ml of aqueous absorbency.
- If the specimen container leaks, the pad is positioned to absorb the spill *inside the specimen bag*.

PEACE OF MIND PERFORMANCE

- When a specimen bag with pre-stuffed absorbent pad is not available, Medegen / Minigrip absorbent pads ARE available – to insert into ANY specimen bag you have.
- Absorbent pads provide the spill containment needed, should the primary specimen container leak inside of the bag.



| Minigrip Code # | Medegen Code # | Description | Qty/cs |
|-----------------|----------------|---|--------|
| ABS3X5MG | ABS3X5MMP | 3" x 5" WHT Super Absorbent Pad 10/100 1000/CS | 1000 |



INTEPLAST GROUP
Inteplast Healthcare

MEDEGEN
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

Customer Service Fax 631-714-6311

E-mail cservice@medegenmed.com