

# EMEGREEN™ EMESIS BAG

#1933

*Convenient personal solution  
that can be easily placed  
anywhere in the facility*

## FEATURES:

**Convenience:** Wide opening is fully supported by a rigid plastic ring to help ensure unobstructed access to the containment bag while providing a circular grasping surface for secure grip and proper positioning around the mouth.

**Capacity:** Deep cylindrical shape holds in excess of 1000 cc/ml (40 fl. oz.) for managing vomit to help prevent spills, and costly clean ups. Durable polyethylene construction is lightweight and tear-resistant.

**Minimize Exposure:** Rigid ring support channels fluid into the containment bag to help reduce spillage and manage exposure.

**Secure Disposal:** Twist, close, and hook bag into plastic ring support notch for easy, one-step disposal of contents and emesis bag. Odors are contained during transport.

**Easy Dispensing:** Wall Hanging Dispenser (#3933D) holds folded containment bags securely and facilitates individual dispensing when needed. Dispensers provided with minimum annual purchase commitment of Emesis Bags, or can be purchased separately. Contact your Medegen Medical Products sales representative for details.

**Areas of Use:** Install virtually anywhere patient services are provided.



**MEDEGEN**  
Medical Products  
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036  
Toll Free 800-511-6298 • Customer Service Fax 631-714-6311  
E-mail [cservice@medegenmed.com](mailto:cservice@medegenmed.com)



1

Always within easy reach when used with the Wall Hanging Dispenser (#3933D).



2

Rigid plastic ring provides a wide opening and secure grasping surface.




3

Twist used Emesis bag to safely contain fluids.



4

Dispose after hooking twisted bag into the support ring notch.



### Emesis Bag (#1933) Specifications

Dimensions:	6.5" X 15" overall, containment bag length 14"
Film:	LDPE, transparent green tint allows content visualization
Graduations:	To 1000cc/ml and 40 fl. oz.
Case Quantity:	100 each/case, 25 bags per sleeve (4)
Case Weight:	3.8 lbs.
Latex Free	

### Wall Hanging Dispenser (#3933D)

Ideal For:  
Emergency Room  
Oncology  
Radiology  
PACU Recovery  
Labor and Delivery  
Pediatrics  
General Medical  
Waiting Areas

Dispenser for use with  
Medegen Medical Products  
part# 1933  
Emesis Bag refills only



For more information, please contact your Medegen Medical Products sales representative at 1-800-511-6298. To view our products online, please visit us at [www.medegenmed.com](http://www.medegenmed.com) or directly with your device by scanning the QR code below.

Medegen Medical Products develops, manufactures and distributes a comprehensive portfolio of high quality disposable medical and surgical products to the healthcare community. Medegen is a leading provider of disposable products, patient bedside products, containment systems for medical waste and protective apparel.



Medegen Medical Products  
209 Medegen Drive  
Gallaway, TN 38036

CUSTOMER SERVICE 800-511-6298  
CUSTOMER SERVICE FAX 631-714-6311  
CUST. SERVICE EMAIL CService@medegenmed.com

INTL CUST. SERVICE IntCS@medegenmed.com  
WEB [www.medegenmed.com](http://www.medegenmed.com)

© 2019 Medegen Medical Products  
LIT# 0066

