

# Surgical & Isolation MASKS

Medegen Medical Products offers a simplified selection of Surgical and Isolation masks. The Surgical masks are ASTM F2100 Level 3 rated and the Isolation/Procedure masks are Level 1 rated.

Our masks are engineered with multiple layers of protection while remaining light, comfortable and providing breathability.

## Surgical masks – 9901 & 9902

- 4 layer of SMS protection
- Fluid Resistant for surgical use
- Ties offer secure fit
- Ear Loops are made with a soft nonirritating material to ensure easy donning and removal

9901, 9902	LEVEL 3
BACTERIAL FILTRATION	98% OR GREATER
SUB-MICRON PARTICULATES FILTRATION EFFICIENT AT .1MICRON	98% OR GREATER
DIFFERENTIAL PRESSURE, mmH <sub>2</sub> O/cm (BREATHABILITY)	5.0 OR GREATER
RESISTANCE TO PENETRATION BY SYNTHETIC BLOOD, MIN PRESSURE IN mmHg FOR PASS RESULTS	160 mmHg
Flame Spread	CLASS 1 (DNI)

## Isolation/Procedure masks

### 9903, 9905 & 9906

- 3 layers of protection
- Ties offer a secure fit
- Ear Loops are made with a soft nonirritating material to ensure easy donning and removal

9903, 9905, 9906	Level 1
BACTERIAL FILTRATION	95% OR GREATER
SUB-MICRON PARTICULATES FILTRATION EFFICIENT AT .1MICRON	95% OR GREATER
DIFFERENTIAL PRESSURE, mmH <sub>2</sub> O/cm (BREATHABILITY)	4.0 OR GREATER
RESISTANCE TO PENETRATION BY SYNTHETIC BLOOD, MIN PRESSURE IN mmHg FOR PASS RESULTS	80 mmHg
Flame Spread	CLASS 1 (DNI)

CODE #	DESCRIPTION	QTY/CS	PACKAGING
9901	Tie-on Surgical Mask, Blue – Level 3	500	50
9902	Ear Loop Surgical Mask, Yellow – Level 3	500	50
9903	Tie-on Procedure Mask, Blue – Level 1	500	50
9905	Ear Loop Procedure Mask, Blue – Level 1	500	50
99906	Fluid Shield Mask, Blue – Level 1	100	25

9901 • Surgical Mask

9902 • Surgical Mask

9905 • Isolation Mask

9903 • Isolation Mask

99906  
Isolation Mask

**MEDEGEN**  
Medical Products  
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

Customer Service Fax 631-714-6311

E-mail [cservice@medegenmed.com](mailto:cservice@medegenmed.com)



**INTEPLAST GROUP**  
Inteplast Healthcare



# Protective APPAREL

## Gowns

Code #	Description	Qty/cs	Packaging
9219	Chemotherapy Gown	25	
99910	Isolation Gown, Yellow	50	10
99911	Impervious Gown, Cuffed, Yellow	50	10
99916	Isolation Gown, Blue, Poly	75	15

## Shoe Covers

Code #	Description	Qty/cs	Packaging
99912	Universal Fit, Non-Conductive, Blue	300	100
99913	Universal Fit, Anti-Skid, Blue	300	100
99914	Extra-Large, Non-Conductive, Blue	300	100

## Masks

Code #	Description	Qty/cs	Packaging
9901	Tie-on Surgical Mask, Blue - Level 3	500	50
9902	Ear Loop Surgical Mask, Yellow - Level 3	500	50
9903	Tie-on Procedure Mask, Blue - Level 1	500	50
9905	Ear Loop Procedure Mask, Blue - Level 1	500	50
99906	Fluid Shield Mask, Blue - Level 1	100	25

## Surgical Caps

Code #	Description	Qty/cs	Packaging
99907	24 in. Bouffant Cap, Blue	500	100
99908	21 in. Bouffant Cap, Blue	500	100
MOBCAP10B	Light Weight 21 in 10gsm Bouffant Cap, Blue	1000	100
MOBCAP14B	Standard Weight 21 in 14gsm Bouffant Cap, Blue	1000	100
99909	Surgeon Cap, Blue	500	100
99922	Surgeon Hood, Blue	500	100



**MEDEGEN**  
Medical Products  
*An Inteplast Group Company*

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

Customer Service Fax 631-714-6311

E-mail [cservice@medegenmed.com](mailto:cservice@medegenmed.com)



**INTEPLAST GROUP**  
*Inteplast Healthcare*