

# Protective HEADWEAR AND SHOE COVERINGS

Medegen Medical Products has you covered from Head to Toe offering various types of personal protective headwear to protective shoe coverings.

Our selection of protective headwear has you covered from non-critical areas to surgical area.

## All of our protective headwear is

- Made from Spunbond polypropylene
- Breathable

## We offer Bouffant Caps that are

- Light weight up to standard weight
- 21" or 24" wide
- Mobcaps are ruched for storage without compromising safety or comfort

## Surgical Caps

- Are breathable at the crown
- Ties to ensure a snug fit

## Surgical Hood


- Full head protection
- Ties at the base of the neck for full protection



99922



99909



MOBCAP

99907



99908



99912

## Bouffant Caps

CODE #	DESCRIPTION	QTY/CS	PACKAGING
99907	24 in. Bouffant Cap, Blue	500	100
99908	21 in. Bouffant Cap, Blue	500	100
MOBCAP10B	Light Weight 21 in. 10gsm Bouffant Cap, Blue	1000	100
MOBCAP14B	Standard Weight 21 in. 14gsm Bouffant Cap, Blue	1000	100

## Surgical Protective Headwear

CODE #	DESCRIPTION	QTY/CS	PACKAGING
99909	Surgeon Cap, Blue	500	100
99922	Surgeon Hood, Blue	500	100

Medegen Medical Products offers shoe covers with and without tread.

Shoe covers offer protection against possible exposure to airborne organisms or contact with contaminated environments throughout a facility.

## Protective Shoe Coverings are made

- From spunbond polypropylene
- With elastic opening to provide a snug fit while providing comfort
- With tread where traction is required

## Protective Shoe Covers

CODE #	DESCRIPTION	QTY/CS	PACKAGING
99912	Universal Fit, Non-Conductive, Blue	300	100
99913	Universal Fit, Anti-Skid, Blue	300	100
99914	Extra-Large, Non-Conductive, Blue	300	100

**MEDEGEN**  
Medical Products  
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

Customer Service Fax 631-714-6311

E-mail [cservice@medegenmed.com](mailto:cservice@medegenmed.com)



**INTEPLAST GROUP**  
Inteplast Healthcare



# Protective APPAREL

## Gowns

Code #	Description	Qty/cs	Packaging
9219	Chemotherapy Gown	25	
99910	Isolation Gown, Yellow	50	10
99911	Impervious Gown, Cuffed, Yellow	50	10
99916	Isolation Gown, Blue, Poly	75	15

## Shoe Covers

Code #	Description	Qty/cs	Packaging
99912	Universal Fit, Non-Conductive, Blue	300	100
99913	Universal Fit, Anti-Skid, Blue	300	100
99914	Extra-Large, Non-Conductive, Blue	300	100

## Masks

Code #	Description	Qty/cs	Packaging
9901	Tie-on Surgical Mask, Blue - Level 3	500	50
9902	Ear Loop Surgical Mask, Yellow - Level 3	500	50
9903	Tie-on Procedure Mask, Blue - Level 1	500	50
9905	Ear Loop Procedure Mask, Blue - Level 1	500	50
99906	Fluid Shield Mask, Blue - Level 1	100	25

## Surgical Caps

Code #	Description	Qty/cs	Packaging
99907	24 in. Bouffant Cap, Blue	500	100
99908	21 in. Bouffant Cap, Blue	500	100
MOBCAP10B	Light Weight 21 in 10gsm Bouffant Cap, Blue	1000	100
MOBCAP14B	Standard Weight 21 in 14gsm Bouffant Cap, Blue	1000	100
99909	Surgeon Cap, Blue	500	100
99922	Surgeon Hood, Blue	500	100



**MEDEGEN**  
Medical Products  
*An Inteplast Group Company*

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

Customer Service Fax 631-714-6311

E-mail [cservice@medegenmed.com](mailto:cservice@medegenmed.com)



**INTEPLAST GROUP**  
*Inteplast Healthcare*