

ISOLATION & PROTECTIVE GOWNS

Gowns offer different levels of protection for various areas of a facility.

Cover Gowns

Made of spunbond polypropylene and are breathable, with ties at the neck and waist with elastic cuffs. Recommended use as a standard gown for Visitors, Food Service, Environmental Service and Housekeeping, where light fluid protection is needed.

Isolation Gowns

Made of polyethylene coated spunbond polypropylene gown, with ties at the neck and waist and has elastic cuffs. Recommended use when there is a low to moderated risk of fluid exposure such as Decontamination, Emergency Room, Food Services or anywhere there is a need for fluid protection.

Fluid Resistant Gown

Made of polyethylene with thumb hook design with apron style neck for easy donning of gown and gloves. Recommend use when there is a greater risk of fluid exposure such as Emergency Room, Decontamination, Food Services or anywhere there is a need for fluid protection.

The benefit of Universal sized gowns is to help facilities minimize stocking multiple gowns in various sizes. Universal size fits most.



99910

CODE #	DESCRIPTION	QTY/CS	PACKAGING
99910	Cover Gown, Universal, Yellow	100	10 /10
99911	Isolation Gown, Universal, Yellow	100	10/10
99916	Fluid Resistant Gown, Universal, Blue	75	5/15



209 Medegen Drive, Gallaway, TN 38036
Toll Free 800-511-6298
Customer Service Fax 631-714-6311
E-mail cservice@medegenmed.com

ACTI-FEND™

Protective APPAREL

Protective Eyewear

Code #	Description	Qty/cs	Packaging
*206-	Heavy Duty Safety Goggles, Clear	12	1
*208-	Safety Glasses, Clear	10	1
*209-	Safety Spectacles, Clear	36	1
*210	Safety Glasses, Clear	12	1
211	Safety Eyeshield, Clear	500	Bulk

*Meets the requirements of ANSI Z87.1

Shoe Covers

Code #	Description	Qty/cs	Packaging
99912	Universal Fit, Non-Conductive, Blue	300	3/100
99913	Universal Fit, Anti-Skid, Blue	300	3/100
99914	Extra-Large, Non-Conductive, Blue	300	3/100

Masks

Code #	Description	Qty/cs	Packaging
9901	Tie-on Surgical Mask, Blue - Level 3	500	10/50
9902	Ear Loop Surgical Mask, Yellow - Level 3	500	10/50
9903	Tie-on Procedure Mask, Blue - Level 1	500	10/50
9905	Ear Loop Procedure Mask, Blue - Level 1	500	10/50
99906	Fluid Shield Mask, Blue	100	4/25

Surgical Caps

Code #	Description	Qty/cs	Packaging
99907	24 in. Bouffant Cap, Blue	500	5/100
99908	21 in. Bouffant Cap, Blue	500	5/100
MOBCAP10B	Light Weight 21 in. 10gsm Bouffant Cap, Blue	1000	10/100
MOBCAP14B	Standard Weight 21 in. 14gsm Bouffant Cap, Blue	1000	10/100
99909	Surgeon Cap, Blue	500	5/100
99922	Surgeon Hood, Blue	500	5/100



MEDEGEN
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

Customer Service Fax 631-714-6311

E-mail cservice@medegenmed.com



INTEPLAST GROUP
Inteplast Healthcare