

---

# PERSONAL **PROTECTION**

---



*Protective Apparel for  
Healthcare Facilities*

# PERSONAL PROTECTION PROD

**Protective apparel protects against possible exposure to airborne organisms, contact with contaminated contaminants throughout a facility from non-critical areas to surgical areas.**

## Protective Headwear

**Bouffant Caps and Bouffant mob caps** are made from spunbond polypropylene and is breathable. We offer a lightweight and a standard weight cap as well as 21" or 24" wide versions. Bouffant mob caps are ruched for compact storage without compromising safety or comfort.

**Surgical head coverings** are made from spunbond polypropylene and are breathable at the crown. Surgical caps ties at the nape of the neck while the Surgical hood offers full head protections and ties at the base of the neck.

### Bouffant Caps

CODE #	DESCRIPTION	QTY/CS	PACKAGING
99907	24 in. Bouffant Cap, Blue	500	5/100
99908	21 in. Bouffant Cap, Blue	500	5/100
MOBCAP10B	Light Weight 21" 10gsm Bouffant Cap, Blue	1000	10/100
MOBCAP14B	Standard Weight 21" 14gsm Bouffant Cap, Blue	1000	10/100

### Surgical Protective Headwear

CODE #	DESCRIPTION	QTY/CS	PACKAGING
99909	Surgeon Cap, Blue	500	5/100
99922	Surgeon Hood, Blue	500	5/100

## Protective Eyewear

We have a variety of protective eyewear in different shapes and designs to meet your specific needs while providing protection from debris, bodily fluids and hazardous items.

- Wrap around eyewear with side shields for added protection
- Clear constructed lenses for unobstructed sight
- Our goggles and glasses fit over most prescription eyewear

CODE #	DESCRIPTION	QTY/CS	PACKAGING
*206-	Heavy Duty Safety Goggles, Clear	12	1
*208-	Safety Glasses, Clear	10	1
*209-	Safety Spectacles, Clear	36	1
*210	Safety Glasses, Clear w/extendable arms	12	1
211	Safety Eye Shield, Clear	500	Bulk

\*Meets the requirements of ANSI Z87.1





# UCTS from Medegen Medical Products

## Isolation and Surgical Masks

Medegen offers a simplified selection of Surgical and Isolation masks. Surgical masks are ASTM F2100 Level 3 rated and the Isolation/Procedure masks are Level 1 rated.

Our masks are engineered with multiple layers of protection while remaining light, comfortable and providing breathability.

### Surgical masks – 9901, 9902

- 4 layers of protection
- Fluid Resistant for Surgical use
- Ties offer secure fit
- Ear Loops are made with a soft non-irritating material to ensure easy donning and removal

### Isolation/Procedure masks – 9903, 9905, 99906

- 3 layers of protection
- Ties offer a secure fit
- Ear Loops are made with a soft nonirritating material to ensure easy donning and removal

9901, 9902	Level 3
Bacterial filtration	98% or greater
Sub-Micron particulates filtration efficient at .1micron	98% or greater
Differential Pressure, mmH <sub>2</sub> O/cm (Breathability)	5.0 or greater
Resistance to penetration by synthetic blood	160 mmHg min
Flame Spread	Class 1 (DNI)
9903, 9905, 99906	Level 1
Bacterial filtration	95% or greater
Sub-Micron particulates filtration efficient at .1micron	95% or greater
Differential Pressure, mmH <sub>2</sub> O/cm (Breathability)	4.0 or greater
Resistance to penetration by synthetic blood	80 mmHg min
Flame Spread	Class 1 (DNI)

CODE #	DESCRIPTION	QTY/CS	PACKAGING
9901	Tie-on Surgical Mask, Blue-Level 3	500	10/50
9902	Ear Loop Surgical Mask, Yellow-Level 3	500	10/50
9903	Tie-on Procedure Mask, Blue-Level 1	500	10/50
9905	Ear-Loop Procedure Mask, Blue-Level 1	500	10/50
99906	Fluid Shield Mask, Blue-Level 1	100	4/25

## Protective Gowns

Gowns offer different levels of protection for various areas of a facility.

**Cover Gowns** – are made of spunbond polypropylene and are breathable, with ties at the neck and waist with elastic cuffs. Recommended use as a standard gown for Visitors, Food Service, Environmental Service and Housekeeping, where light fluid protection is needed.

**Isolation Gowns** – are made of polyethylene-coated spunbond polypropylene gown, with ties at the neck and waist and has elastic cuffs. Recommended use when there is a low to moderated risk of fluid exposure such as Decontamination, Emergency Room, Food Services or anywhere there is a need for fluid protection.

**Fluid Resistant Gown** – Made of polyethylene with thumb hook design with apron style neck for easy donning of gown and gloves. Recommend use when there is a greater risk of fluid exposure such as Emergency Room, Decontamination, Food Services or anywhere there is a need for fluid protection.

*The benefit of Universal sized gowns is to help facilities minimize stocking multiple gowns in various sizes. Universal size fits most.*


CODE #	DESCRIPTION	QTY/CS	PACKAGING
99910	Cover Gown, Universal, Yellow	100	10/10
99911	Isolation Gown, Universal, Yellow	100	10/10
99916	Fluid Resistant Gown, Universal, Blue	75	5/15

## Protective Shoe Covers

Shoe covers are made of spunbond polypropylene with elastics opening providing a snug fit with comfort. Our Anti-Skid shoe covers offer tread where traction is required. Non-Conductive shoe covers are available to control static electricity.

CODE #	DESCRIPTION	QTY/CS	PACKAGING
99912	Universal Fit, Non-Conductive, Blue	300	3/100
99913	Universal Fit, Anti-Skid, Blue	300	3/100
99914	Extra-Large, Non-Conductive, Blue	300	3/100

# SEE ALL THE GREAT MEDEGEN OFFERINGS IN OUR PRODUCT CATALOGS



To view our catalogs  
online, please visit us at  
[www.medegenmed.com](http://www.medegenmed.com)  
or directly with your  
device by scanning the  
QR code below.



360 Motor Pkwy, Suite 800 • Hauppauge, NY 11788 • 800.511.6298  
Customer Service Fax: 631.714.6311 • Customer Service: [CSservice@medegenmed.com](mailto:CSservice@medegenmed.com)

[www.medegenmed.com](http://www.medegenmed.com)

LIT# 0097

