



HOSPITAL FURNITURE

*Quality, Comfort and
Convenience*



PATIENT BEDSIDE CABINETS

Our patient Bedside Cabinets come assembled and are offered in wood grain laminate or ABS plastic.

- Color offered in wood laminate: Maple
- Color offered in ABS: Blue
- ABS plastic Bedside Cabinet top has a lip on 3 sides preventing item from falling off
- Left or Right mountable hanger for either Bedside Cabinet
- Wood laminate comes with a slide out shelf, draw, closet and a vented drawer
- ABS plastic comes with a slide out shelf, drawer, closet and a vented shelf
- 4 x 2" casters, front casters are locking
- Dimension: L 22" x W 19" x H 30"



PATIENT OVERBED TABLE

Medegen's Overbed Tables are designed to be:

- STABLE – with an “H” shaped base and four casters (2 x 1" and 2 x 2")
- DURABLE – we offer two options: wood grain laminate or ABS Plastic
- FUNCTIONAL – our Overbed Table has a gas spring lifting column with a fingertip self-locking handle enabling Healthcare workers to adjust the table height from 30" up to 41". Accommodating the Patient Bed, Stretcher, Chair/Recliner or a Wheelchair
- Overbed Table Dimensions are: L 35" x W 17"
- Overbed Tables are available in ABS plastic color Blue and Wood Laminate color Light Maple
- Table Weight capacity is 60 LBS

HOSPITAL FURNITURE

RECLINER-SLEEPER

5 Operating Handles

Operating Handles are on the left and right side of the recliner, making it easy to move or assisting with converting the recliner to a sleeper.



2 Safety Knobs

Safety Knobs are located on the left and right sides of the recliner-sleeper. Either knob can easily convert the recliner to a sleeper and back to a recliner.



1

2

3



4

5

6

4 Gas Spring Back

The recliner has a Gas Spring back making the transition from recliner to sleeping surface effortless.



1 Arm Rests

Arm rests are made of soft comfortable foam.



6 Braking Caster

Recliner comes with 4 x 3" casters, 2 are locking and are toe operated.



3 Sleeper

Recliner-Sleeper has an ergonomic sleeping surface with drop down armrests providing an unobstructed sleeping area.



Medegen Recliner-Sleepers are offered in 3 different sizes and 2 colors:

SMALL

Chair: L 34" x W 24" x H 39"
Bed: L 77" x W 24" x H 18"

Sleeping surface: L 73" x W 21"

MEDIUM

Chair: L 34" x W 27" x H 39"
Bed: L 79" x W 27" x H 18"

Sleeping surface: L 75" x W 22"

LARGE

Chair: L 35" x W 30" x H 42"
Bed: L 85" x W 30" x H 18"

Sleeping surface: L 80" x W 24"

BROWN


BEIGE

RECLINER-SLEEPER *ECHL Series*



1 Safety Knobs

Safety Knobs are located on the left and right sides of the recliner-sleeper. Either knob can easily convert the recliner to a sleeper and back to a recliner.



2 Operating Handles

Operating Handles are on the left and right side of the recliner, making it easy to move or assisting with converting the recliner to a sleeper.



4 Sleeper

The recliner has a mechanical spring making the transition from recliner to sleeping surface effortlessly.

5 Sleeping Surface

Recliner-Sleeper has an ergonomic sleeping surface with a built-in head rest.



6 Casters


Recliners come with 6 casters, 4 x 1.6" at the base of the recliner and 2 x 3" casters at the rear of the recliner for easy transportation of recliner.

3 Storage

Recliner Sleeper has storage in the seat of the recliner. Storage has handles for easy opening and has hinges for safe opening and closing.

7 Arm Rests

Arm rests are stationary and are covered in PVC with soft comfortable foam.



Medegen Recliner-Sleeper series is offered in 3 different colors:

Recliner-Sleeper – ECHL Series overall dimensions

Recliner Length: 37.8" Width: 30.7" Height: 37.8"

Overall Sleeper Dimensions: Length: 76" Width: 23.6"

Sleeping Surface Length: 74.8" Width: 23.67"

BROWN

BLACK

BEIGE

To Order: 800.511.6298

www.medegenmed.com

RECLINER-SLEEPER *ECGM Series*



1 Safety Knobs

Safety Knobs are located on the left and right sides of the recliner-sleeper. Either knob can easily convert the recliner to a sleeper and back to a recliner.



2 Operating Handles

Operating Handles are on the left and right side of the recliner, making it easy to move or assisting with converting the recliner to a sleeper.



3 Storage


Recliner Sleeper has storage in the seat of the recliner. Storage has handles for easy opening and has hinges for safe opening and closing.

4 Sleeper

The recliner has a mechanical spring making the transition from recliner to sleeping surface effortlessly.

5 Sleeping Surface

Recliner-Sleeper has an ergonomic sleeping surface with drop down armrests providing an unobstructed sleeping area.



6 Arm Rests

Arm rests are made of soft comfortable foam.



7 Casters

Recliners come with 6 casters, 4 x 1.6" at the base of the recliner and 2 x 3" casters at the rear of the recliner for easy transportation of recliner.

Medegen Recliner-Sleeper series is offered in 3 different colors:

Recliner-Sleeper – ECGM Series overall dimensions

Recliner Length: 36.41" Width: 27.17" Height: 36.6"

Sleeping Surface Length: 74.8" Width: 21.7"

BROWN

BLACK

BEIGE



FDA



OTHER GREAT PRODUCTS FROM MEDEGEN


the **KALIKART**™ series of **HEALTHCARE CARTS**



See our full line of carts in this informative trifold product brochure!



CHEVALIER Hydraulic STRETCHER



See all the features in our comprehensive trifold product brochure!



To view our products online, please visit us at www.medegebenmed.com or directly with your device by scanning the QR code to the right.



360 Motor Pkwy, Suite 800 • Hauppauge, NY 11788 • 800.511.6298
Customer Service Fax: 631.714.6311 • Customer Service: CSservice@medegenmed.com
International Customer Service: IntCSservice@medegenmed.com

www.medegebenmed.com

LIT# 0092

