

Your Source For
POST MORTEM
BAGS
and Kits

Post Mortem Bags are an integral part of the health / care continuum and are designed for the respectful and safe transport of deceased individuals.

Medegen® brand post mortem bags (PMBs) are thoughtfully designed to ensure ease of use and containment assurance.

Features & Benefits

- **Expansive Opening:** Wide “D” shape or perimeter zipper opening for simplified patient loading and closure.
- **Rugged Zipper Closure:** Heavy-duty, zipper closure with extra-long zipper pulls for easy grasping and fast closure.
- **Contrasting Color:** White bag with thick black zipper for fast open / closure identification.
- **Leak Resistant:** Wide-welded side seams and liquid tight zipper simplifies containment safety.
- **Identification:** Featuring an exclusive external document pouch with hook & loop closure and identification tags for the range of identification needs.
- **Writable Surface:** For convenience; also printed with Universal Precautions Statement.
- **Materials:** Offered in durable DEHP-Free PVC and Chlorine-Free PEVA versions.
- **Complete Size Range:** From infant (18" x 28") to larger individuals (48" x 92"), we have you covered.



PMB3690D

MEDEGEN
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298 • Customer Service Fax 631-714-6311

E-mail cservice@medegenmed.com



To view our products online, please visit us at www.medegenmed.com or directly with your device by scanning the QR code to the left.

Your Source For **Post Mortem Bags & Kits**

Code #	Size (WxL) (in)	Size	Material	Zipper Access Style	Film Color	ID Kit*	Post Mortem Kit**	Qty/cs
PMB3690D	36x90	Adult	PVC Embossed	D-Style	White	✓		10
PMB3690DK	36x90	Adult	PVC Embossed	D-Style	White	✓	✓	10
PMB3690P	36x90	Adult	PVC Embossed	Perimeter	White	✓		10
PMB3690DPE	36x90	Adult	PEVA Embossed	D-Style	White	✓		10
PMB3690DPEK	36x90	Adult	PEVA Embossed	D-Style	White	✓	✓	10
PMB4892D	48x92	LARGE Adult	PVC Embossed	D-Style	White	✓		10
PMB2848D	28x48	Child	PVC Embossed	D-Style	White	✓		20
PMB1828D	18x28	Infant	PVC Embossed	D-Style	White	✓		30

***ID Kit:** (3) ID tags with string attached to each, (1) 4" black zip tie packed in zip closure bag for personal effects.

****Post Mortem Kit:** (1) poly bag for personal items (12" x 8" x 20"), (2) 60" ties, (3) 36" ties, (1) nonwoven chin strap, (2) underpads 16.5" x 23.5".

