

Collection & **CONTAINMENT**



Today, more than ever, medical professionals can't afford to settle for just average performance in their health care products. It's critical for you to know you can count on every product performing as expected. That's why Medegen Medical Products pays meticulous attention to every detail of its full line of Collection & Containment Solutions. Working under the most stringent quality controls ensures that every product from Medegen Medical Products is crafted to exacting standards, from design and manufacture to packaging and delivery. Superior performance begins with uncompromising quality. Let Medegen help you provide...

Better Healthcare Through Better Products.



Table of Contents

Patient Belonging Bags	2-6
Lab Specimen Transport Bags	7-11
Zip Closure Bags	12-16
Specialty Medical Bags	17-18
Sterility Maintenance Covers	19
Medical Equipment and Cart Covers	20-21
Liners	21-37
Bedside Product Liners	21
Biohazardous Waste and Chemo Waste Bags	23-26
Autoclavable Biohazardous Waste Bags	27-28
Laundry and Linen Bags	29-31
Water Souble Laundry Bags	32
Hampers and Accessories	33
Open Top Bags	34
Institutional Trash Can Liners	35-37
Sharps Containers	38-42
Stackable Sharps Containers	39-40
Non-Stackable Sharps Containers	41
Sharps Hardware	42
Spill Control Solutions	43-45
Solidifiers	43
Spill Clean-up Kits	44-45
Protective Apparel	46
Gowns	46
Protective Eyewear	46
Index	47-49

Patient Belonging, Lab Specimen and Zip Bags



Medegen Medical Products is your BEST source for cost saving, innovative patient care, medical specimen, and healthcare packaging needs.

- **Patient Belonging Bags** – Check out our Industry leading line of patient belonging and setup bags, including the economical TIE&GO® Bag!
- **Lab Specimen Bags** – Now with industry leading Lab Guard® and Tear Zone® brands Medegen offers the broadest range of easy-to-use, secure specimen transport packaging.
- **Zip Bags** – Because not everything is healthcare or the lab has a custom package, Medegen offers a broad range of zip bag sizes with thick, secure white zippers for easy isolation and protective packaging.

Our field Sales Consultants understand the unique demands of health care providers and can help in selecting the appropriate products for intended applications, all with the goal to help you optimize your resources and operations efficiencies.

Patient Belonging & Setup Bags

Patient Belongings Bags are a convenient and cost effective means to collect, transport and store patient belongings from the time of admission, transfer to room, between departments and ultimate discharge. Bags come in a variety of sizes, styles and strengths to handle most personal belongings.

- Choice of handles/closures include: Drawstring, Draw Tape, Rigid Handle with snap lock feature, TIE&GO handle, Zip Closure and Adhesive Tabs
- Durable polyethylene construction with built in gusset expands capacity, is tear resistant and provides a stable platform for adding/removing clothing
- Colored film for personal privacy/security or clear for loss prevention
- Standard labeling facilitates an easy to recognize write-on space for client name and room number. Also available without print
- Ideal solution to help organize and manage patient belongings or set-up materials during stay
- Set-up bags identify specific use applications and client specific ease of control
- **Customized Printing** provides immediate identification and public relations/brand awareness for your facility (name, logo, address, etc.) or specific department (ER, Same Day Surgery, Maternity, etc.). *Minimum order quantity required*
- Typical areas of use: Admissions, ER, Outpatient Clinics, Labor and Delivery, Discharge and wherever client clothing/personal effects are collected



Patient Belongings Bags — LLDPE Film, Flat Pack, Drawstring

Code#	Bag/Print Color	Size (inches)	Capacity	Usage	Qty/cs	Packaging
Print/Label: Patient Belongings, Patient Name, Room Number.						
DSPB20C	Clear/Blue	20 X 18 X 3	Standard	Regular Duty	250	25
DSPB20W	White/Blue	20 X 18 X 3	Standard	Regular Duty	250	25
3539	White/Blue	20 X 19 X 4	Large	Heavy Duty	250	25
3537	Clear/White & Blue	20 X 20 X 3	Large	Heavy Duty	250	25
50-20	White/Blue	20 X 20 X 4	Large	Regular Duty	250	25
50-50	Clear/White & Blue	20 X 20 X 4	Large	Regular Duty	250	25
50-53	Green/White	20 X 20 X 4	Large	Regular Duty	250	25
50-80	Fuchsia/Black	20 X 24 X 4	Extra Large	Heavy Duty	250	25
DSPB2020B	Blue/Black	20 X 20 X 4	Large	Regular Duty	250	25
3549	White/Blue	28 X 32 X 4	Extra Large	Heavy Duty	150	15



DSPB20W




DSPB2020B



50-50



50-53



50-80

Patient Belongings Bags

★ Indicates Made in the USA



Patient Belongings Bags — LLDPE Film, Flat Pack, Rigid Handle

Code#	Bag/Print Color	Size (inches)	Capacity	Usage	Qty/cs	Packaging
Print/Label: Patient Belongings, Patient Name, Room Number.						
RHPB20C	Clear/Blue	20 x 18 x 3	Standard	Regular Duty	250	25
RHPB20W	White/Blue	20 x 18 x 3	Standard	Regular Duty	250	25
50-21	White/Blue	20 x 20 x 4	Large	Regular Duty	250	25
70-21	White/Blue	20 x 20 x 4	Large	Heavy Duty	250	25
50-51	Clear/White & Blue	20 x 20 x 4	Large	Regular Duty	250	25
70-51	Clear/White & Blue	20 x 20 x 4	Large	Heavy Duty	250	25
70-31	Blue/Black	22 x 19 x 3	Extra Large	Extra Heavy Duty	250	25



Patient Belongings Bags - LLDPE Film, Flat Pack, Drawtape

Code#	Bag/Print Color	Size (inches)	Capacity	Usage	Qty/cs	Packaging
Print/Label: Personal Belongings, Patient Room, Room Number.						
DTPB20W	White/Blue	20 x 18 x 3	Standard	Regular Duty	250	25



TIE&GO® Patient Belonging Bag with Loop Handles

Code#	Bag/Print Color	Size (inches)	Capacity	Usage	Qty/cs	Packaging
Print/Label: Personal Belongings, Patient Room, Room Number.						
TBPB20W ★	White/Blue	16 x 8 x 26	Large	Regular Duty	500	50

Patient Set-Up Bags with Draw Tape

Code#	Bag/Print Color	Size (inches)	Usage	Qty/cs	Packaging
Print/Label: Patient Name, Room Number, Time/Date Issued, Equipment in Use, Hospital Property Do Not Remove (English/Spanish).					
50-30	Clear/White & Black	12 x 16	Regular Duty	500	50



Patient Set-Up Bags with Drawstring

Code#	Bag/Print Color	Size (inches)	Usage	Qty/cs	Packaging
Print/Label: Patient Name, Room Number, Time/Date Issued, Equipment in Use, Hospital Property Do Not Remove (English/Spanish)					
UFRSU1216	Clear/White & Black	12 x 16	Regular Duty	500	50

Patient Belongings Bags

★ Indicates **Made in the USA**

Traditional T-Shirt Bag - Loop Handles

Minigrip Code#	Bag/Print Color	Size (inches)	Capacity	Usage	Qty/cs
Print/Label: THANK YOU					
THW1A ★	White/Red	12 x 3 x 21	Small		1000



THW1A

Patient Garment Bag with Zipper and Hanger

Minigrip Code#	Bag/Print Color	Size (inches)	Capacity	Usage	Qty/cs
Print/Label: None					
VGB-HPL	White/None	24 x 40			50



VGB-HPL

Patient Care and Utility Bags - LLDPE Film, Flat Pack, Die Cut Handle

Code#	Bag/Print Color	Size (inches)	Capacity	Usage	Qty/cs
Print/Label: None					
SG9WHI	White/None	9 x 12	Small	Regular Duty	1000
SG15WHI	White/None	15x18+4	Small	Regular Duty	500
SG18WHI	White/None	18x18+4	Medium	Regular Duty	500

Patient Care and Utility Bags - LLDPE Film, Flat Pack, Reinforced Foldover Die Cut Handle

Code#	Bag/Print Color	Size (inches)	Capacity	Usage	Qty/cs
Print/Label:★ THANK YOU					Δ
FO9WHI	White/None	9 X 12	Small		1000
FO12WHI	White/None	12 x 15	Small		500
FO12THA*	White/Red	12 x 15	Small		500


Patient Specialty Care Bags



50-45

Patient Bedside with Adhesive Tabs (HDPE)

Code#	Bag/Print Color	Size (inches)	Qty/cs	Packaging
Print/Label: Blue Floral Print				
50-45	White/Blue	7 x 3.5 x 11.5	2000	100



3660

Patient Assistive Device Bag - Zip Closure

Code#	Bag/Print Color	Size (inches)	Qty/cs	Packaging
Print/Label: Patients Hearing Aids, Dentures, Glasses				
3660	Fuchsia/Clear & Black	7.5 x 13	1000	100



UFPVE912

Patient Valuables Bag Tamper Evident Adhesive Closure

Code#	Bag/Print Color	Size (inches)	Qty/cs	Packaging
Print/Label: Sequential Numbered Receipt, Writable Patient/Content Block				
UFPVE912	Opaque Gray/Black	9 x 12	1000	100

Lab Specimen Transport Bags

Laboratory Specimen Transport Bags facilitate collection, storage, identification, and transport of laboratory specimens and cultures within the hospital or surgery center, as well as to off-site facilities and laboratories. Bags contain a dual pouch system which separates and protects samples from the accompanying paperwork.

- **Lab Guard®**  **Branding**, the specimen packaging trusted by hospitals and labs for over 25 years. Also featuring popular Spec-Zip® and SAFEGUARD® brands
- **TearZone® EZ-Open Feature** makes specimen bags faster to open than the zipper when a specimen is in processing
- **Expanded Temperature / Stat Selector Check Boxes** for even greater versatility when using Medegen's Lab Guard specimen bags
- Available in a variety of standard sizes in clear or colored film with easy to use zip closure, reclosable adhesive, or permanent tamper evident adhesive closure systems
- Typical areas of use: Phlebotomy Labs, Chemistry & Hematology, Pathology, Toxicology, Microbiology, Emergency Room, Treatment Areas, Drug Screening, or wherever Routine Clinical Laboratory Specimen Collections are drawn
- **FDA Registered!**



4914



SBL2X69B5



4510

Lab Guard® **Speci-Zip® & SAFEGUARD® Lab Specimen Transport Bags — Zip Closure**

Lab Guard®, Spec-Zip® & SAFEGUARD® brands provide convenient, easy-to-use, cost-saving features you just can't get with any other bag in a range of combinations. Choose below from single zip and liquid tight double zip closures, standard or cost saving destroyable biohazard symbol, TearZone® EZ open feature or NO TearZone, paperwork pouch flap or NO flap, STAT+Temp or just Temp selector boxes and more!

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SBL2X46B		4x6	Single Zip	Clear/Blk & Org		X	Standard		no flap	DispBag	20/50	1000
IP46B3T		4x6	Single Zip	Clear/Blk & Org		X	Standard	Score	w/flap	Zip Bag	10/100	1000
	3670	4x6	Single Zip	Clear/Blk & Org	X		Standard		w/flap	Zip Bag	10/100	1000
	49-91	6x6	Single Zip	Clear/Blk & Red	X		Standard		w/flap	Zip Bag	10/100	1000
SBL2X69B5		6x9	Double Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
SBL2X69B		6x9	Double Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBOX	2/500	1000
IP69B		6x9	Single Zip	Clear/Blk & Red			Standard	Score	w/flap	Zip Bag	10/100	1000
	4075	6x9	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4914	6x9	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4914.ORG	6x9	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	no flap	Zip Bag	10/100	1000
	D4078	6x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4510	7x9	Slider Zip	Clear/Blk & Org	X		Standard		w/flap	Zip Bag	10/50	500
SBL2X88B		8x8	Single Zip	Clear/Blk & Org		X	Standard		no flap	DispBag	20/50	1000
	49-95	8x8	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X810B5		8x10	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
SBL2X810B		8x10	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBOX	2/500	1000
	4915	8x10	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4915.ORG	8x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X912B		9x12	Single Zip	Clear/Blk & Org		X	Standard	TearZone®	no flap	DispBag	20/50	1000
	49-99	9x12	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4925	10x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X118B		11x18	Single Zip	Clear/Blk & Org		X	Standard		no flap	Zip Bag	20/50	1000
SBL2X1215B		12x15	Single Zip	Clear/Blk & Org		X	Standard		no flap	Zip Bag	40/25	1000
MGTZ1215		12x15	Single Zip	Clear/Blk & Org		X	Standard	TearZone®	no flap	DispBag	5/100	500
IP1215B3T		12x15	Single Zip	Clear/Blk & Org		X	Standard	Score	w/flap	Zip Bag		250
	4916	12x15	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	5/100	500
	59-94	12x13x4	Single Zip	Clear/Blk & Red	NONE	NONE	Standard		w/flap	Zip Bag	10/100	1000
ZLABTRANSPT		12.75x12x6	Single Zip	Clear/Blk & Org	NONE	NONE	Standard		no doc pch	Zip Bag	5/100	500
ZLAB1420		14x20	Single Zip	Clear/Blk & Org	NONE	NONE	Standard		no doc pch	Zip Bag	5/100	500
IP2024B3T		20x24	Single Zip	Clear/Blk & Org		X	Standard		w/flap	Zip Bag	1/100	100

Lab Specimen Transport Bags

Lab Guard® & Speci-Zip® Lab Specimen Transport Bags — Zip Closure, STAT

Lab Guard® & Speci-Zip STAT Specimen Bags are available in yellow, red, and blue color schemes. Other available features include cost-saving TearZone®, destroyable biohazard symbol or standard biohazard symbol and DispenserBag® or standard zip bag packaging.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
	3785	6x9	Single Zip	Yel Opaque/Blk		X	Standard	TearZone®	w/flap	Zip Bag	10/100	1000
MGTZ69S		6x9	Single Zip	Red Tint/Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
IP69RSTAT2		6x9	Single Zip	Red Tint/Blk		X	Standard	Score	w/flap	Zip Bag	10/100	1000
	4001	6x9	Single Zip	Red Tint/Blk		X	Standard	TearZone®	w/flap	Zip Bag	10/100	1000
IP69RSTAT		6x9	Single Zip	Red Opaque/Blk		X	Standard	Score	w/flap	Zip Bag	10/100	1000
SBL2X69S		6x9	Single Zip	Clr/Red,Wht,Blk		X	Standard		no flap	DispBag	20/50	1000
	4002	6x9	Single Zip	Blu Opaque/Blk		X	Standard	TearZone®	w/flap	Zip Bag	10/100	1000
IP810RSTAT		8x10	Single Zip	Red Opaque/Blk		X	Standard	Score	w/flap	Zip Bag	10/100	1000



MGTZ69S



SBL2X69S



4001



3785



4002



SBL2AP69B

Lab Guard® & Speci-Zip® Lab Specimen Transport Bags — Zip Closure with Absorbent Pad Insert (3x5 Pad)

These Lab Guard® & Speci-Zip Specimen Bags feature a super-absorbent 3x5 absorbent pad to soak up leaking specimen cup fluids.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SBL2AP69B		6x9	Double Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
IP69BAS		6x9	Single Zip	Clear/Blk & Red			Standard	Score	w/flap	Zip Bag	10/100	1000
SBL2AP810B		8x10	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
SBL2AP1215B		12x15	Single Zip	Clear/Blk & Org		X	Standard		no flap	DispBag	20/50	1000

Lab Guard® Lab Specimen Transport Bags — Zip Closure, Utility (No Biohazard Print)

Lab Guard® UTILITY Specimen Bags are not printed with a biohazard symbol for hospitals and labs with in-facility processing and no transportation outside of the facility.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
	49-90	6x6	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000
SBL2X69		6x9	Double Zip	Clear/Blk & Org	NONE	NONE	NONE	TearZone®	no flap	DispBag	20/50	1000
	49-96	6x9	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000
	49-93	8x8	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000

49-90



SBL2X69



Lab Specimen Transport Bags

Lab Guard® & Speci-Zip® Lab Specimen Transport Bags — Zip Closure, Color-Codes

Lab Guard® & Speci-Zip® Color-Code Specimen Bags feature a range of color blocks on the front of bag to help hospitals and labs segregate specimens with specific processing needs. Features cost-saving TearZone®, destroyable biohazard symbol.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
MGTZ69BLU	4811B	6x9	Single Zip	Clear/Blu,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
IP69BBL		6x9	Single Zip	Blue/Red,Blk			Standard	Score	w/flap	Zip Bag	10/100	1000
MGTZ69YEL		6x9	Single Zip	Clear/Yel,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
IP69BY		6x9	Single Zip	Yellow/Red,Blk			Standard	Score	w/flap	Zip Bag	10/100	1000
MGTZ69GRN	4811G	6x9	Single Zip	Clear/Grn,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
IP69BG		6x9	Single Zip	Green/Red,Blk			Standard	Score	w/flap	Zip Bag	10/100	1000
MGTZ69RD		6x9	Single Zip	Clear/Red,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
MGTZ69PUR		6x9	Single Zip	Clear/Pur,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000



MGTZ69RED



MGTZ69BLU



MGTZ69YEL



MGTZ69GRN



MGTZ69PUR



Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure

Speci-Gard® Specimen Bags feature the innovative, easy-to-use, liquid-tight PRESS AND CLOSE® permanent adhesive closure, TearZone® EZ open system, DispenserBag® convenience packaging, and in some cases, a cost-saving destroyable biohazard symbol.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
95-CMS59		5x9	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-66SPBIO		6x6	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	no pouch	DispBag	100	2000
UF95-600BIO		6x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-600BPM		6x10	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	DispBag	100	1000
95-CMS810		8x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-610BIO		10x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-1318BIOA		13x18	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	Zip Bag	50	500
UF95-1517BIO		15x17	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	Zip Bag	50	500



UF95-600BIO



UF95-600BPM



UF95-600STB



UF95-600STBR

To Order: 800.511.6298 • www.medegenmed.com

Lab Specimen Transport Bags

Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure - STAT

Speci-Gard® STAT Specimen Bags are available in red and yellow film colors with all featuring the easy-to-use, liquid-tight, PRESS AND CLOSE® permanent adhesive closure, TearZone® EZ open system, cost saving destroyable biohazard symbol and DispenserBag® convenience packaging.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SPEC-PNYRT		6x10	Adhesive	Red Tint/Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
UF95-600STBR		6x10	Adhesive	Red Opaque/Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
UF95-600STB		6x10	Adhesive	Yel Opaque/Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000

Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure & Absorbent Pad Insert (3x5 Pad)

Speci-Gard® Specimen Bags with absorbent pads are available with either single adhesive pouch (UF95-BDES) or DUAL-adhesive pouch (SPECI-CC). Item SPECI-CC secures both the specimen and the paperwork in separate pouches so that neither can be tampered with, without destroying the bag (chain-of-custody).

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
UF95-BDES		6x10	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	DispBag	10/100	1000
SPECI-CC		6x10	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	Adhesive	DispBag	10/100	1000



UF95-BDES



SPEC-PNYB



SPEC-PNYG



SPEC-PNYP

Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure - Color Codes


Speci-Gard® Color-Code PRESS AND CLOSE® Adhesive Closure Specimen Bags feature a range of film colors to help hospitals and labs segregate specimens with specific processing needs.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SPEC-PNYB		6x10	Adhesive	Blu Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
UF95-600BBL		6x10	Adhesive	Blu Opq/Blk & Red		X	Standard	TearZone®	w/flap	DispBag	10/100	1000
SPEC-PNYG		6x10	Adhesive	Gm Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
SPEC-PNYP		6x10	Adhesive	Pur Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
SPEC-PNYY		6x10	Adhesive	Yel Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000

SPEC-PNYY



UF95-600BBL



UF95-600



Lab Specimen Transport Bags

Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure, Utility (No Biohazard Print)

Speci-Gard® UTILITY Specimen Bags are not printed with a biohazard symbol for hospitals and labs with in-facility processing and no transportation outside of the facility. All bags offer TearZone® EZ-Open Feature and DispenserBag® Convenience Packaging.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
UF95-66SP		6x6	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	10/100	2000
UF95-600		6x10	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	10/100	1000
UF95-610		10x10	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	10/100	1000

Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure, UV Protection

Speci-Gard® style PRESS AND CLOSE® Permanent Adhesive Closure Bag features opaque gray white film that prevents light from entering the bag, protecting light sensitive contents. Access to contents is available with EZ-Tear Perforation.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SPECI-UV		6x10	Adhesive	Wht Opq/Blk & Red	NONE	NONE	Standard	Perforation	w/flap	DispBag	10/100	1000

Lab Specimen Transport Bags: Permanent Adhesive Closure, Foldover Closure

Lab specimen bags feature a traditional foldover permanent adhesive closure and score line for opening during processing.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
	59-89	6x9	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	10/100	1000
	4060	6x10	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	10/100	1000
	59-93	10x10	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	10/100	1000

Lab Specimen Transport Bags: Adhesive Closure, Specialty Applications

Item 1130 features rugged welded side seams and permanent adhesive closure for pneumatic tube systems.

Item 39-97 features a reclosable adhesive closure for worry-free post-closure access.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
	1130	6x9	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Tear Notch	no flap	Bag	10/100	1000
	39-97	6x10	Reclosable Adh	Clear/Blk & Red	NONE	NONE	Standard	Pull Open	no flap	Zip Bag	10/100	1000



SPECI-UV



4060



1130



39-97

Minigrip® Red Line™ Zip Closure Bags

- Featuring full gauge thickness and the thickest, most rugged zippers on the market.
- Sizes 10" in length and above feature liquid-tight double zipper for enhanced content containment.



Minigrip® Red Line™ Zip Closure Bags - Clear, LDPE Film, Flat Pack

Minigrip Code #	Bag/Print Color	Size (inches)	Packaging	Qty/cs	Qty/Master cs
2 Mil Plain					
MGRL2P0202	Clear/None	2 x 2	100	1,000	10,000
MGRL2P0203	Clear/None	2 x 3	100	1,000	10,000
MGRL2P0205	Clear/None	2 x 5	100	1,000	5,000
MGRL2P0208	Clear/None	2 x 8	100	1,000	10,000
MGRL2P0303	Clear/None	3 x 3	100	1,000	5,000
MGRL2P0304	Clear/None	3 x 4	100	1,000	5,000
MGRL2P0305	Clear/None	3 x 5	100	1,000	5,000
MGRL2P0306	Clear/None	3 x 6	100	1,000	
MGRL2P0404	Clear/None	4 x 4	100	1,000	5,000
MGRL2P0406	Clear/None	4 x 6	100	1,000	5,000
MGRL2P0407	Clear/None	4 x 7	100	1,000	5,000
MGRL2P0408	Clear/None	4 x 8	100	1,000	
MGRL2P0505	Clear/None	5 x 5	100	1,000	
MGRL2P0507	Clear/None	5 x 7	100	1,000	
MGRL2P0508	Clear/None	5 x 8	100	1,000	
MGRL2P0510*	Clear/None	5 x 10	100	1,000	
MGRL2P0604	Clear/None	6 x 4	100	1,000	
MGRL2P0606	Clear/None	6 x 6	100	1,000	
MGRL2P0608	Clear/None	6 x 8	100	1,000	
MGRL2P0609	Clear/None	6 x 9	100	1,000	
MGRL2P0610*	Clear/None	6 x 10	100	1,000	
MGRL2P0612*	Clear/None	6 x 12	100	1,000	
MGRL2P0808	Clear/None	8 x 8	100	1,000	
MGRL2P0810*	Clear/None	8 x 10	100	1,000	
MGRL2P0812*	Clear/None	8 x 12	100	1,000	
MGRL2P0912*	Clear/None	9 x 12	100	1,000	
MGRL2P1012*	Clear/None	10 x 12	100	1,000	
MGRL2P1013*	Clear/None	10 x 13	100	1,000	
MGRL2P1212*	Clear/None	12 x 12	100	1,000	
MGRL2P1215*	Clear/None	12 x 15	100	1,000	
MGRL2P1315*	Clear/None	13 x 15	100	1,000	
MGRL2P1318*	Clear/None	13 x 18	100	1,000	

2 Mil Writable White Block

MGRL2W0203	Clear/None	2 x 3	100	1,000	10,000
MGRL2W0304	Clear/None	3 x 4	100	1,000	5,000
MGRL2W0305	Clear/None	3 x 5	100	1,000	5,000
MGRL2W0406	Clear/None	4 x 6	100	1,000	5,000
MGRL2W0508	Clear/None	5 x 8	100	1,000	
MGRL2W0609	Clear/None	6 x 9	100	1,000	
MGRL2W0810*	Clear/None	8 x 10	100	1,000	
MGRL2W0912*	Clear/None	9 x 12	100	1,000	


4mil Plain with Hang Hole

MGRL4PH02503	Clear/None	2.5 x 3	100	1,000	
MGRL4PH02506	Clear/None	2.5 x 6	100	1,000	
MGRL4PH0304	Clear/None	3 x 4	100	1,000	
MGRL4PH0404	Clear/None	4 x 4	100	1,000	
MGRL4PH0406	Clear/None	4 x 6	100	1,000	
MGRL4PH0408	Clear/None	4 x 8	100	1,000	
MGRL4PH0508	Clear/None	5 x 8	100	1,000	
MGRL4PH0604	Clear/None	6 x 4	100	1,000	
MGRL4PH0606	Clear/None	6 x 6	100	1,000	
MGRL4PH0608	Clear/None	6 x 8	100	1,000	
MGRL4PH1012*	Clear/None	10 x 12	100	500	

*These items come with double zippers



DOUBLE ZIPPER



SINGLE ZIPPER



MGRL2P0609

GreenLine™ Zip Closure Bags

Minigrip® Red Line® Zip Closure Bags - Clear, LDPE Film, Flat Pack

(continued)

Minigrip Code #	Bag/Print Color	Size (inches)	Packaging	Qty/cs
4mil Plain No Print				
MGRL4P0609	Clear/None	6 x 9	100	1,000
MGRL4P0612*	Clear/None	6 x 12	100	1,000
MGRL4P0808	Clear/None	8 x 8	100	1,000
MGRL4P0810*	Clear/None	8 x 10	100	1,000
MGRL4P0912*	Clear/None	9 x 12	100	1,000
MGRL4P1010*	Clear/None	10 x 10	100	500
MGRL4P1013*	Clear/None	10 x 13	100	500
MGRL4P1212*	Clear/None	12 x 12	100	500
MGRL4P1215*	Clear/None	12 x 15	100	500
MGRL4P1318*	Clear/None	13 x 18	100	500
MGRL4P1320*	Clear/None	13 x 20	50	250
MGRL4P1424*	Clear/None	14 x 24	50	250
4mil Writable White Block with Hang Hole				
MGRL4WH0304	Clear/None	3 x 4	100	1,000
MGRL4WH0305	Clear/None	3 x 5	100	1,000
MGRL4WH0406	Clear/None	4 x 6	100	1,000
MGRL4WH0508	Clear/None	5 x 8	100	1,000
MGRL4WH0609	Clear/None	6 x 9	100	1,000
MGRL4WH0810*	Clear/None	8 x 10	100	1,000
MGRL4WH0912*	Clear/None	9 x 12	100	1,000
6mil Plain No Print				
MGRL6P0406	Clear/None	4 x 6	50	500
MGRL6P0608	Clear/None	6 x 8	50	500
MGRL6P0810*	Clear/None	8 x 10	50	500
MGRL6P0912*	Clear/None	9 x 12	50	500
MGRL6P1012*	Clear/None	10 x 12	50	500
MGRL6P1215*	Clear/None	12 x 15	50	250
MGRL6P1318*	Clear/None	13 x 18	50	250


*These items come with double zippers



MGRL4PH0404



MGRL4WH0406



MGRL4P1010

Minigrip® GreenLine™ - Biodegradable Zip Closure Bags Zip Closure Bags - Clear, LDPE Film, Flat Pack

- The same performance as traditional Minigrip® Red Line™ bags and tested to biodegrade 15% within 1 year per ASTM D5511 test results.
- Sizes 10" in length and above feature liquid-tight double zipper for enhanced content containment.
- Biodegradable - Every part of the bag can biodegrade - film, zipper material - EVERYTHING!
- Proven ASTM D5511 tested by an independent lab to be biodegradable.

Minigrip Code #	Bag/Print Color	Size (inches)	Packaging	Qty/cs	Qty/Master cs
2 Mil Plain, No Print					
MGBD2P0202	Clear/None	2 x 2	100	1,000	10,000
MGBD2P0203	Clear/None	2 x 3	100	1,000	10,000
MGBD2P0304	Clear/None	3 x 4	100	1,000	5,000
MGBD2P0406	Clear/None	4 x 6	100	1,000	5,000
MGBD2P0507	Clear/None	5 x 7	100	1,000	
MGBD2P0609	Clear/None	6 x 9	100	1,000	
MGBD2P0912*	Clear/None	9 x 12	100	1,000	

*These items come with double zippers

MGBD2P0912



MGBD2P0406

To Order: 800.511.6298 • www.medegenmed.com

Zip Closure Bags



Zip Closure Bags provide a convenient low cost solution for kitting, for collection, storage and transport of personal items or pre-packaging of individual articles for client use.

- **Robust White Zipper** is easy to open and truly effortless to shut, keeping your product securely inside
- **Dispenser Bag® Convenience Packing** provides neat storage and serves one bag at a time. Available for most bag sizes
- Available in a wide range of sizes with and without writable white block to provide multi-use versatility throughout your facility
- Typical areas of use: Pharmacy, Central Supply, Kitting, Dietary and throughout the facility

Reloc-Zippit® Zip Closure Bags – Clear, LDPE Film, Flat Pack

The premier economy zip bag brand opens easily and features a popular thick white zipper for fast, secure closure and content containment.




Minigrip Code #	Medegen Code#	Bag/Print Color	Size (inches)	Qty/cs	Packaging
2mil Plain No Print					
MGZ2P0202	Z2.0202	Clear /None	2 x 2	1,000	100
MGZ2P0203	Z2.0203	Clear /None	2 x 3	1,000	100
MGZ2P0208	-	Clear /None	2 x 8	1,000	100
MGZ2P0303	Z2.0303	Clear /None	3 x 3	1,000	100
MGZ2P0304	Z2.0304	Clear /None	3 x 4	1,000	100
MGZ2P0305	Z2.0305	Clear /None	3 x 5	1,000	100
MGZ2P0306	-	Clear /None	3 x 6	1,000	100
MGZ2P0404	Z2.0404	Clear /None	4 x 4	1,000	100
MGZ2P0406	Z2.0406	Clear /None	4 x 6	1,000	100
MGZ2P0407	-	Clear /None	4 x 7	1,000	100
MGZ2P0408	Z2.0408	Clear /None	4 x 8	1,000	100
MGZ2P0505	-	Clear /None	5 x 5	1,000	100
MGZ2P0507	Z2.0507	Clear /None	5 x 7	1,000	100
MGZ2P0508	Z2.0508	Clear /None	5 x 8	1,000	100
MGZ2P0510	-	Clear /None	5 x 10	1,000	100
MGZ2P0604	-	Clear /None	6 x 4	1,000	100
MGZ2P0606	Z2.0606	Clear /None	6 x 6	1,000	100
MGZ2P0608	Z2.0608	Clear /None	6 x 8	1,000	100
MGZ2P0609	Z2.0609	Clear /None	6 x 9	1,000	100
MGZ2P0610	Z2.0610	Clear /None	6 x 10	1,000	100
MGZ2P0708	-	Clear /None	7 x 8	1,000	100
MGZ2P0808	Z2.0808	Clear /None	8 x 8	1,000	100
MGZ2P0810	Z2.0810	Clear /None	8 x 10	1,000	100
MGZ2P0812	-	Clear /None	8 x 12	1,000	100
MGZ2P0912	Z2.0912	Clear /None	9 x 12	1,000	100
MGZ2P1010	-	Clear /None	10 x 10	1,000	100
MGZ2P1012	Z2.1012	Clear /None	10 x 12	1,000	100
MGZ2P1013	Z2.1013	Clear /None	10 x 13	1,000	100
MGZ2P1212	Z2.1212	Clear /None	12 x 12	1,000	100
MGZ2P1215	Z2.1215	Clear /None	12 x 15	1,000	100
MGZ2P1315	Z2.1315	Clear /None	13 x 15	1,000	100

Zip Closure Bags

Reloc-Zippit® Zip Closure Bags – Clear, LDPE Film, Flat Pack (continued)

Minigrip Code #	Medegen Code#	Bag/Print Color	Size (inches)	Qty/cs	Packaging
MGZ2P1318	Z2.1318	Clear /None	13 x 18	1,000	100
MGZ2P1420	-	Clear /None	14 x 20	500	100
MGZ2P1824	-	Clear /None	18 x 24	250	25
MGZ2P2020	-	Clear /None	20 x 20	500	50
2 Writable White Block					
MGZ2W0203	Z2.0203W	Clear /None	2 x 3	1,000	100
MGZ2W0304	Z2.0304W	Clear /None	3 x 4	1,000	100
MGZ2W0305	Z2.0305W	Clear /None	3 x 5	1,000	100
MGZ2W0406	Z2.0406W	Clear /None	4 x 6	1,000	100
MGZ2W0508	Z2.0508W	Clear /None	5 x 8	1,000	100
MGZ2W0609	Z2.0609W	Clear /None	6 x 9	1,000	100
MGZ2W0810	-	Clear /None	8 x 10	1,000	100
MGZ2W0912	Z2.0912W	Clear /None	9 x 12	1,000	100
MGZ2W1215	-	Clear /None	12 x 15	1,000	100
4mil Plain No Print					
MGZ4P02503	-	Clear /None	2.5 x 3	1,000	100
MGZ4P0304	Z4.0304	Clear /None	3 x 4	1,000	100
MGZ4P0305	-	Clear /None	3 x 5	1,000	100
MGZ4P0404	-	Clear /None	4 x 4	1,000	100
MGZ4P0406	Z4.0406	Clear /None	4 x 6	1,000	100
MGZ4P0508	Z4.0508	Clear /None	5 x 8	1,000	100
MGZ4P0606	-	Clear /None	6 x 6	1,000	100
MGZ4P0608	Z4.0608	Clear /None	6 x 8	1,000	100
MGZ4P0609	Z4.0609	Clear /None	6 x 9	1,000	100
MGZ4P0808	Z4.0808	Clear /None	8 x 8	1,000	100
MGZ4P0810	Z4.0810	Clear /None	8 x 10	1,000	100
MGZ4P0912	Z4.0912	Clear /None	9 x 12	500	100
MGZ4P1010	-	Clear /None	10 x 10	1,000	100
MGZ4P1012	Z4.1012	Clear /None	10 x 12	500	100
MGZ4P1212	Z4.1212	Clear /None	12 x 12	500	100
MGZ4P1215	-	Clear /None	12 x 15	500	100
	Z4.1215	Clear /None	12 x 15	1000	100
MGZ4P1318	Z4.1318	Clear /None	13 x 18	500	100
MGZ4P1424	-	Clear /None	14 x 24	250	50
MGZ4P1824	-	Clear /None	18 x 24	250	25
MGZ4P2020	-	Clear /None	20 x 20	250	50
4mil Writable White Block					
MGZ4W0305	-	Clear/None	3 x 5	1000	100
MGZ4W0406	-	Clear/None	4 x 6	1000	100
MGZ4W0508	-	Clear/None	5 x 8	1000	100
MGZ4W0609	-	Clear/None	6 x 9	1000	100
MGZ4W0912	-	Clear/None	9 x 12	500	100
6mil Plain No Print					
MGZ6P1318	-	Clear /None	13 x 18	250	50

Food Service



COLORZIP® Food Service Double Zip Closure

Food Service professionals rely on Minigrip recloseable bags with **COLORZIP™** double zipper. From sandwich to jumbo 2 gallon size, you'll never run short of options when it comes to your food storage needs. Signature **blue tint bags** are designed to be tougher, stronger and more liquid tight than ever - perfect for storage and short term freezing. Includes write-on white block area for easy content & date identification. Packed in sturdy dispenser box packaging.



MGDZW175G

Minigrip Code #	Bag/Print Color	Size (inches)	Capacity	Qty/cs
MGSZW115S50*	Clear/None	6.5 x 6	Sandwich	500
MGDZW175Q	Blue Tint/None	7 x 8	Quart Storage	500
MGDZW270Q	Blue Tint/None	7 x 8	Quart Freezer	300
MGDZW175G	Blue Tint/None	10.5 x 11	Gallon Storage	250
MGDZW270G	Blue Tint/None	10.5 x 11	Gallon Freezer	200
MGDZW175J	Blue Tint/None	13 x 15.5	2 Gal Storage	100
MGDZW270J	Blue Tint/None	13 x 15.5	2 Gal Freezer	100

* Single Zip Closure



Slider Zip Closure Bags

Enjoy the convenience of the easy slide open/close feature that's made the Minigrip® **slider zip** bag **SO POPULAR!** These seal-top bags are built to tackle a wide variety of applications, especially in environments where gloves are worn. They open and close quickly and easily the first time everytime, reducing unnecessary repetitive motion. Packed in sturdy **dispenser box** packaging*.



MGS270J

Minigrip Code #	Bag/Print Color	Size (inches)	Capacity	Packaging*	Qty/Master CS
MGS270Q	Clear/None	8 x 7	Quart Freezer	250	1000
MGS270L	Clear/None	9 x 12	Letter	250	1000
MGS270G	Clear/None	10.5 x 11	Gallon Freezer	250	1000
MGS270J	Clear/None	14.5 x 16	2 Gal Freezer	250	500

Specialty Medical Bags

UV Protection Zip Bags (Amber Color)

Take comfort in knowing the integrity of critical light sensitive products won't be compromised during testing, shipping and storage. Contents can be viewed through the film. Film contains UV blockers to prevent bag contents from harmful light impact. Premium thick, white zipper opens quickly and easily. Complies with USP procedure 66-1 with light transmission not exceeding 10% at any wavelength. No Print.

Minigrip Code#	Bag Color	Size (inches)	Packaging	Qty/CS	Qty/Master CS
MGUV3P0203	Amber	2 x 3	100	1,000	10,000
MGUV3P02509	Amber	2.5 x 9	100	1,000	4,000
MGUV3P0305	Amber	3 x 5	100	1,000	4,000
MGUV3P0406	Amber	4 x 6	100	1,000	2,000
MGUV3P0508	Amber	5 x 8	100	1,000	2,000
MGUV3P0608	Amber	6 x 8	100	1,000	
MGUV3P0808	Amber	8 x 8	100	1,000	
MGUV3P0814	Amber	8 x 14	100	500	
MGUV3P1212	Amber	12 x 12	100	500	



MGUV3P0814

UV Protection Adhesive Closure Bag (White/Gray Opaque)

Speci-Gard® style Press and Close® Adhesive Closure Bag features opaque gray white film that prevents light from entering the bag, protecting light sensitive contents. Access to contents is available with EZ-Tear Perforation. See page 11, item SPECI-UV.

Anti-Static Zip Bags (Pink Tint Color)


- Bags qualify as Type 2 under MIL-B-81705D entitled barrier material flexible and electrostatic free
- Bags are tested to dissipate an induced 5000 volt charge to technical zero in accordance with FTS 101C, method 4046.1
- Smaller sizes packaged in DispenserBag® convenience packaging
- Pink color film indicates anti-static properties; thick, durable white zip closure. No Print.

Minigrip Code#	Bag Color	Size (inches)	Qty/cs	Packaging
MGST4P02503	Pink Tint	2.5 x 3	1,000	100
MGST4P0305	Pink Tint	3 x 5	1,000	100
MGST4P0406	Pink Tint	4 x 6	1,000	100
MGST4P0608	Pink Tint	6 x 8	1,000	100
MGST4P0609	Pink Tint	6 x 9	1,000	100
MGST4P0810	Pink Tint	8 x 10	1,000	100
MGST4P0912	Pink Tint	9 x 12	500	50
MGST4P1212	Pink Tint	12 x 12	500	50
MGST4P1318	Pink Tint	13 x 18	500	50



MGST4P1212

Specialty Medical Bags



UF95-BL10

Disposal Bag For Used IV Blood Bags

Secure permanent adhesive closure bags designed specifically for disposal of used IV blood bags. Instructions for use and a biohazard logo are printed on the outside of the pocket. Two pockets for blood bag and paperwork.

Minigrip Code#	Bag Color	Size (inches)	Qty/cs	Packaging
UF95-BL10	Clear	10 x 10	1000	100



MGCHEM069

Chemotherapy Transport Bags - LDPE Film, Flat Pack, Zip Closure

For safe transport of chemotherapy drugs. These secondary storage containers are used by doctors, hospitals and pharmacies to transport chemo medications from pharmacy to patient and return of empty containers to pharmacy.

Minigrip Code#	Medegen Code#	Bag/Print Color	Size	Additional	Qty/cs	Packaging
MGCHEM069		Clear/Yel & Blk	6x9	English/Spanish	1000	100
MGCHEM0912		Clear/Yel & Blk	9x12	English/Spanish	500	100
MGCHEM01215		Clear/Yel & Blk	12x15	English/Spanish	500	100
	9218	Yellow/Black	12x15	English/Spanish;CYTA	250	50
	47-CBS1	Yellow/Black	12x15	English/Spanish;CYTA	100	



47-CBS1



UFEMB1010

Emesis Containment Bag - LDPE Film, Flat Pack, Permanent Adhesive Closure

Wide mouth and side handles make this bag very user-friendly. OPAQUE white film hides contents. Once used, liquid-tight adhesive closure provides optimal infection and odor protection. Biohazard print for infectious contents.

Minigrip Code#	Medegen Code#	Bag/Print Color	Size	Qty/cs	Packaging
UFEMB1010		White/Org & Blk	10x10	1000	100

Sterility Maintenance Covers

Sterility Maintenance Covers protect sterile packs from dust, dirt, odors, moisture and foreign materials, and are used to package, store, transport and maintain a dust free environment. **Instrument Covers and Scope Transport Bags** facilitate transport and easy identification of contaminated and clean devices. Typical areas of use: Central Processing and anywhere medical instruments and sterile products are packaged/processed.

Sterility Maintenance Covers — LDPE Film, Flat Pack, Self-Seal Adhesive Strip

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Qty/cs	Packaging
Print/Label: Tear and Seal Instructions.					
810	Clear/Blue	6 x 10	2.25	1000	250
815	Clear/Blue	8 x 12	2.25	1000	250
820	Clear/Blue	10 x 15	2.25	1000	250
825	Clear/Blue	12 x 20	2.25	500	250
830	Clear/Blue	16 x 22	2.25	500	250
835	Clear/Blue	16 x 30	2.25	250	—
840	Clear/Blue	24 x 30	2.25	250	—
845	Clear/Blue	25 x 37	2.25	100	—



- Self Seal Adhesive strip secures contents and provides a tamper-evident barrier for sterile integrity
- Easy opening linear tear strip allows for quick access and eliminates the risk associated with using scissors or blades to open

Instrument Transport Covers — LDPE Film, Flat Pack, Twist Tie Closures

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Qty/cs	Packaging
Print/Label: Red Biohazard, Contaminated Materials, Instrument Transport Bag.					
52-04	Clear/Red	15 x 5 x 31	1.30	250	25



Scope Transport Bags - LDPE Film, Flat Pack, Zip Closure with Die Cut Handle

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Qty/cs
Print/Label: Scope Transport Bag, Contaminated /Dirty.				
5050C	Red/Black	24 x 20	4.0	200
Print/Label: Scope Transport Bag, Processed/Clean.				
5051	Clear/Black	24 x 20	4.0	200



PROTEX-RAY™ – X-RAY Cassette Covers

X-Ray and Imaging cassettes come in contact with a high number of patients and can be a source of infection for both patients and staff. Use of disposable single use Cassette Covers prevents infection transmission. Generous sizing and foldover flap closure.

Minigrip Code#	Bag Color	Size (inches)	Qty/cs	Packaging
UFXCC1012	Clear	15 x 12.5	500	100
UFXCC1417	Clear	20 x 17.5	500	100



UFXCC1012

Medical Equipment and Cart Dust Covers

★ Indicates **Made in the USA**



Equipment Dust Covers are designed for non sterile protection of medical equipment, linen, and supplies during transport or while not in use. Joint Commission recommends the use of protective covers for shielding clean linen from dirt and contaminants.

- Effective dust and dirt barrier helps keep contents cleaner longer
- Covers are puncture resistant, easy to position, and come in a wide range of sizes
- Available in flat pack or rolls for easy dispensing
- Typical areas of use: Central Processing and anywhere medical instruments/equipment are used

Medical Equipment and Cart Dust Covers — LLDPE Film, Coreless Roll

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Qty/cs	Qty/roll
Print/Label: Equipment Cover.						
43-01	★	Clear/Black	16 x 14 x 38	1.35	150	15
43-05	★	Clear/Black	28 x 22 x 60	1.35	50	10
Print/Label: Biohazard Symbol with Contaminated /Dirty.						
X2740	★	Red/Black	38 x 26 x 48	.90	50	10

Medical Equipment and Cart Dust Covers — LLDPE Film, Cored Roll

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Qty/cs	Qty/roll
No Print/Label.						
6220	★	Clear/None	26 x 22 x 60	1.50	50	50
6490	★	Clear/None	36 x 24 x 43	1.0	300	300/Disp. Box
48-07	★	Clear/None	38 x 26 x 48	.90	50	50/Disp. Box
48-10	★	Clear/None	38 x 26 x 48	1.85	50	50
T1133	★	Clear/None	38 x 26 x 60	1.50	50	50
48-19	★	Clear/None	45 x 27 x 48	.90	50	50/Disp. Box
48-20	★	Clear/None	45 x 27 x 48	1.85	50	50
48-29	★	Clear/None	55 x 27 x 62	.90	50	50
48-30	★	Clear/None	55 x 27 x 62	1.85	50	50
6625	★	Yellow/None	55 x 27 x 62	1.0	50	50
D6624	★	Clear/None	55 x 27 x 62	2.0	50	50
48-39	★	Clear/None	64 x 27 x 67	.90	50	50
48-40	★	Clear/None	64 x 27 x 67	1.85	50	50
6092	★	Clear/None	64 x 27 x 67	2.0	50	50

Medical Equipment and Cart Dust Covers

★ Indicates **Made in the USA**

Medical Equipment and Cart Dust Covers - LLDPE Film, Flat Packed

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Qty/cs
Print/Label: Equipment Cover				
6089	★ Clear/Black	7 x 8	4.0	1,000
6058	★ Clear/Black	16 x 22	2.0	1,000
6059	★ Clear/Black	24 x 30	2.0	500
43-03	★ Clear/Black	16 x 14 x 56	1.35	100
6456	★ Gray/Black	38 x 27 x 50	1.30	100

Bed Covers - LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Qty/cs
No Print/Label				
5139	Clear/None	75 x 125	1.0	100

Bedside Product Liners

As a leading manufacturer of bedpans, wash basins and other patient bedside essentials, Medegen introduces **BEDSIDE PRODUCT LINERS**.

The liners are designed to be a **time-saving & hygienic tool** for nurses and caregivers. It takes time to rinse and sanitize bedpans, commode buckets, wash basins and other items once they've been used. Medegen bedside product liners simply install into the bedside product and once used, are easily removed, tied shut and discarded. And the bedside product isn't soiled! Bedside product liners offer **peace-of-mind performance**, giving caregivers more time to spend with patients and not cleaning devices.

Bedpan & Commode Liner - LLDPE Film, Flat Seal, Drawtape Closure, Absorbent Pad Insert, Roll Pack

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Qty/cs	Packaging
Print: No Print.					
BPL2214DT-W	White	22 x 14	1.40	200	20



Liners



Medegen liners and bags are manufactured in the USA to meet specified performance standards in accordance with Federal and specific state guidelines. Medegen produces a broad range of liners and bags to meet your specific applications needs:

LLDPE and HDPE Can Liners Kick Bucket Liners
Biohazard Red Waste Bags Laundry & Linen Bags
Biohazard Autoclavable Waste Bags Water Soluble Linen Bags
 Hamper Stands

- Premium LLDPE, HDPE, PP and PVA resins to ensure consistent, superior strength performance
- Latex free and environmentally safe for incineration disposal
- Stocked in multiple colors to facilitate waste stream and laundry segregation. Custom prints and sizes also available. Contact us for a quote!,
- **Biohazard** Bags are available that comply with Federal interstate and specific state requirements for 165 gram dart impact resistance (ASTM D1709) and 480 grams tear strength (ASTM D1922) and US DOT CFR 49 173-197.

Almost 100% of our liners are made in the USA for your peace of mind.

Look for the ★ next to items made in the USA.



MEDEGEN
 Medical Products
 An Intoplast Group Company

To Order: 800.511.6298 • www.medegenmed.com

Biohazardous Waste Bags

Biohazardous Waste Collection Bags are used to collect, store, and transport biohazardous / infectious waste within the health care facility, outpatient surgery center, alternate care clinic, physician's office and lab.

- Offered in both Linear Low Density (LLDPE) and High Density (HMW HDPE) resins. Bags < 2 mil are star-sealed while those 2 mil and greater are flat sealed. Compact coreless perforated or interleaved rolls provide compact storage and ease of use.
- Bags can be CUSTOM printed with facility or specific department name. Contact us for a quote!
- Biohazardous / Infectious Waste labeling in both English and Spanish for easy identification of hazardous waste to reduce the risk of contamination during collection and transport.
- Meets 165 gram of ASTM D1709 Dart Impact Resistance.
- Dual-tested bags* meet the required 165 gram ASTM D1709 Dart Impact Resistance and the 480 gram ASTM D1922 Elmendorf Tear Strength properties. Dual tested bags are required for use in Florida, Pennsylvania, and California and for Federal, US DOT interstate transport.
- Accessory hampers may be used with biohazard bags for convenient placement throughout the facility.

★ Indicates **Made in the USA**

Biohazardous Waste Bags - Meet A.S.T.M. Dart Impact/Elmendorf Tear & DOT CFR 49 173-197 Requirements*, LLDPE Film, Star Seal, Coreless Rolls

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.						
50-40	Red/Black	11 x 14	1.50	1-3	500	100**
50-42	Red/Black	11 x 14	1.50	1-3	500	50**
2634 ★	Red/Black	19 x 23	1.25	5-7	400	20
F116 ★	Red/Black	23 x 23	1.20	7-10	500	25
37-97 ★	Red/Black	24 x 24	1.20	8-10	150	150
45-50 ★	Red/Black	24 x 24	1.25	8-10	250	25
F124 ★	Red/Black	24 x 30	1.20	8-12	500	25
45-54 ★	Red/Black	24 x 32	1.25	12-16	250	25
2635 ★	Red/Black	25 x 35	1.35	12-16	200	20
45-30 ★	Red/Black	28 x 31	1.25	20	250	25
F145 ★	Red/Black	30 x 36	1.20	20-30	250	25
F117 ★	Red/Black	31 x 41	1.20	20-30	250	25
47-71 ★	Red/Black	31 x 41	1.60	20-30	100	10
F118 ★	Red/Black	33 x 39	1.20	33	250	25
F130 ★	Red/Black	37 x 50	1.30	44	150	15
F110 ★	Red/Black	37 x 50	1.50	44	150	15
F109 ★	Red/Black	38 x 45	1.20	44	200	20
CF2359 ★	Red/Black	38 x 45	1.75	44	125	25
F149 ★	Red/Black	40 x 45	1.20	40-45	250	25
F135 ★	Red/Black	40 x 46	1.20	40-45	100	10
45-43 ★	Red/Black	40 x 46	1.25	40-45	150	15
45-45 ★	Red/Black	40 x 55	1.25	45-55	150	15

**Flat seal, flat pack



To Order: 800.511.6298 • www.medegenmed.com

MEDEGEN
Medical Products
An Inteplast Group Company

Biohazardous Waste Bags

★ Indicates Made in the USA

Biohazardous Waste Bags - Meet A.S.T.M. Dart Impact/Elmendorf Tear & DOT CFR 49 173-197 Requirements*, LLDPE Film, Flat Seal, Coreless Rolls

Code#Bag/Print Color			Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol							
RD610	★	Red/Black	6 x 5 x 14.25	1.50	1-3	200	20
RD630	★	Red/Black	7.5 x 7 x 19	1.25	2-3	200	20
RD650	★	Red/Black	8.5 x 7.5 x 24	1.25	5-6	200	20
47-41	★	Red/Black	23 x 8 x 41	2.40	20-30	100	10
47-43	★	Red/Black	23 x 17 x 46	2.40	40-45	100	10
F169	★	Red/Black	38 x 45	3.0	40-45	100	10
47-45	★	Red/Black	23 x 17 x 55	2.40	45-55	100	10
45-75	★	Red/Black	55 x 60	2.0	60	50	10

Chemotherapy Waste Bags - CYTA Symbol, Meet A.S.T.M. Dart Impact/Elmendorf Tear & DOT CFR 49 173-197 Requirements, LLDPE Film, Flat Seal Coreless Rolls (2 & 3mil), Flat Packed (4mil)

Chemotherapy Waste Collection Bags are used for transport of anti neoplastic medications and chemotherapy waste. Bags feature reinforced seals for leak resistance and meet labeling requirements for chemotherapy waste collection.

- Color Coded and Bold Labeling with International Biohazard Symbol for easy identification
- Heavy Duty Construction for Puncture Resistance
- For Chemo Drug Transport Bags, See Page 18

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Chemotherapy Waste, CYTA Symbol.						
9233	 Yellow/Black	24 x 34	4.0	12-16	100	flat pack
9230	 Yellow/Black	25 x 34	2.0	12-16	250	25
4210	 Yellow/Black	30 x 43	3.0	20-30	100	10
47-CBL4	 Yellow/Black	31 x 41	4.0	20-30	100	flat pack
X3009B	 Yellow/Black	45 x 49	1.5	50	100	10

Biohazardous Waste Bags - Meet A.S.T.M. Dart Impact Test Requirements, Zip Closure, Flat Pack

Code#	Bag/Print Color	Size (inches)	Gauge/Mil	Capacity/gal	Qty/cs	Packaging
Print/Label: English/Spanish Biohazard, Biohazard Symbol						
4923	Red/Black	9 x 12	2.0	0.75	1000	100
4930	Red/Black	12 x 12	2.0	1	1000	100
4921	Red/Black	12 x 15	2.0	1.5	500	100

Dressing Disposal Bag

Biohazard dressing disposal bag is packed 50 to a dispenser clip. All bags are perforated for easy removal from the clip. Clip packs can be hung in the room, on a cart or kept neatly in a drawer.

Minigrip Code#	Bag Color	Size (inches)	Qty/cs	Packaging
UFDD1114	Red	11 x 14	500	50



4923
(print)




UFDD1114

Biohazardous Waste Bags

Biohazardous Waste Bags - Meet A.S.T.M. Dart Impact Test Requirements, Star Seal, Coreless Rolls (<2mil)

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol.							
115H	★	Red/Black	23 X 23	1.10	7-10	500	25
116	★	Red/Black	23 x 23	1.50	7-10	500	25
57-00	★	Red/Black	24 x 24	1.05	8-10	500	25
108M	★	Red/Black	25 x 34	1.20	12-16	250	25
D2210	★	Red/Black	30 x 36	1.0	20-30	250	25
117M	★	Red/Black	30.5 x 41	1.20	30	250	25
117	★	Red/Black	30.5 x 41	1.50	30	250	25
2020	★	Red/Black	31 x 38	1.85	30	100	10
47-55	★	Red/Black	30.5 x 41	1.10	30	250	25
104M	★	Red/Black	33 x 40	1.20	33	250	25
2045	★	Red/Black	38 x 45	1.0	44	100	10
172M	★	Red/Black	38 x 45	1.25	44	250	25
109	★	Red/Black	38 x 45	1.50	44	200	20
110	★	Red/Black	38 x 45	1.75	44	125	25
D2209	★	Red/Black	40 x 46	1.20	40-45	250	25
47-73	★	Red/Black	40 x 46	1.60	40-45	100	10
2122	★	Red/Black	43 x 55	1.50	55	200	20



Biohazardous Waste Bags - Meet A.S.T.M. Dart Impact Test Requirements, Flat Seal, Coreless Rolls (2mil and greater)

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol							
2302	★	Red/Black	14 X 19	2.0	2-3	500	25
2305	★	Red/Black	25 x 35	2.25	12-16	200	20
2202	★	Red/Black	30 x 40	3.0	30	100	10
117A	★	Red/Black	30.5 x 41	2.0	30	250	25
119	★	Red/Black	30.5 x 41	3.0	30	125	25
2204	★	Red/Black	38 x 45	2.0	44	100	10
169	★	Red/Black	38 x 45	3.0	44	100	10
173	★	Red/Black	43 x 48	2.0	50	100	10
2125	★	Red/Black	43 x 55	3.0	55	100	20

Biohazardous Waste Bags

★ Indicates Made in the USA



Anti-Static Kick Bucket Liner - NFPA 56A, Twist Ties, LLDPE Film, Star Seal, Coreless Rolls

Anti-Static Kick Bucket Liners comply with Operating Room requirements of NFPA Standard 56A. Bags contain a carbon additive to help prevent the buildup of electrostatic charges. 50-02 for round buckets (fits Medegen 58130 and 58160 and other 12" rounds); 50-04 for oval buckets.

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
50-02	★	Red/Black	12 x 8 x 24	1.30	7	500	25
50-04	★	Red/Black	18 x 8 x 16	1.30	7	500	50




Biomedical Waste Containers — Corrugated Box with Liner, Flat Pack

Biomedical Waste Containers provide an effective and inexpensive means for handling large amounts of biomedical waste while controlling disposal costs. Marked UN-3291.

- 3-Part System: Outer Corrugated Container, Inner Corrugated Container, Plastic Liner separation
- Locking Lid, Puncture and Leak Resistant, Incinerator-Safe. Store flat to minimize space

Code#		Bag/Print Color	Size (inches)	Capacity/gal	Qty/cs
10-2011	★	Red/Black & Red	12 x 24 x 12	11gal/18.9L	14



44-05
(print)



RS404816RH
(print)

Biohazardous Waste Bags - Meet A.S.T.M. Dart Test Requirements, HDPE Film, Star Seal, Coreless Roll

Code#		Bag/Print Color	Size (inches)	Gauge/Mic	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Biohazardous/Infectious Waste, Biohazard Symbol. (RS codes English, Spanish, French)							
44-03		Red/Black	17 x 18	8.00	4-6	1,000	100
RS242411RH	★	Red/Black	24 x 24	11.00	8-10	1,000	50
RS243311RH	★	Red/Black	24 x 33	11.00	12-16	1,000	50
RS304314RH	★	Red/Black	30 x 43	14.00	30-32	250	25
44-05	★	Red/Black	31 x 41	16.00	30-32	250	25
HRD334014	★	Red/Black	33 x 40	14.00	33	500	25
44-01	★	Red/Black	33 x 40	16.00	33	250	25
RS384616RH	★	Red/Black	38 x 46	16.00	44	250	25
44-13	★	Red/Black	40 x 46	16.00	40-45	200	20
RS404816RH	★	Red/Black	40 x 48	16.00	40-45	250	25

Autoclavable Biohazardous Waste Bags

Autoclavable Biohazardous Waste Bags are built to withstand the rigors of Steam Sterilization while meeting the requirements for the collection, storage, transport, and processing of biohazardous waste within the acute care facility. Autoclavable bags made in the USA from high strength Polypropylene (PP) Film have an operating temperature of 274°F (134°C) and can tolerate temperatures of up to 285°F (140°C). **Dual-Tested Autoclave Bags** are **ASTM D1922** Elmendorf tear tested to exceed 480g of tear strength (MD, TD) and **ASTM D1709** tested to exceed 165g of dart impact. Standard autoclave bags are tested to exceed 165g of dart impact. Bags contain an integral steam processing indicator that turns brown when autoclaved.

- Extreme temperature resistance and the BEST available tear and dart impact strength make the Dual-Tested Autoclave Bags your "peace of mind" choice. Complies with US DOT 49CFR Sections 173 – 197 and regulations in all 50 states including those with dual-test bag requirements.
- Standard autoclave bags are tested to exceed 165g of dart impact for compliance in the majority of US states (without dual-test bag requirements).
- Biohazardous or Infectious labeling in both English and Spanish or French
- International Biohazard Symbol in easy to see contrasting color
- Available in several color-coded styles for easy identification by department
- Twist ties included for easy secure closure
- Custom prints and sizes too. Contact us for a quote!

 Indicates **Made in the USA**

Dual Tested Autoclavable Biohazardous Waste Bags with Indicator, Meet A.S.T.M. Dart Impact/Elmendorf Tear & DOT CFR 49 173-197 Requirements - Flat Seal, Coreless Roll, Twist Ties

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Max. Temp 285F Print/Label: Biohazard Labeling, Biohazard Symbol, English/Spanish.						
AC1419C 	Clear/Black	14 x 19	2.0	3-4	500	25
AC1419R 	Red/Black	14 x 19	2.0	3-4	500	25
AC1923C 	Clear/Black	19 x 23	2.0	7-10	200	20
AC1923R 	Red/Black	19 x 23	2.0	7-10	200	20
AC2530C 	Clear/Black	25 x 30	2.0	12-14	200	20
AC2530R 	Red/Black	25 x 30	2.0	12-14	200	20
AC2535C 	Clear/Black	25 x 35	2.0	13-16	200	20
AC2535R 	Red/Black	25 x 35	2.0	13-16	200	20
AC3138C 	Clear/Black	31 x 38	2.0	20-30	100	10
AC3138R 	Red/Black	31 x 38	2.0	20-30	100	10
AC3847C 	Clear/Black	38 x 47	2.0	40-50	100	10
AC3847R 	Red/Black	38 x 47	2.0	40-50	100	10




Autoclavable Biohazardous Waste Bags - Polypropylene Film with Indicator, Meet A.S.T.M. Dart Impact Requirements - Flat Seal, Coreless Roll w/Twist Ties

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Max. Temp 285°F. Print/Label: Biohazard Labeling, Biohazard Symbol, English/Spanish.						
880 	Orange/Black	8 x 12	1.80	1-2	400	20
5001.4 	Red/Black	8 x 12	2.0	1-2	400	20
852 	Clear/Red	12 x 24	1.80	2-3	200	20
5010 	Red/Black	12 x 24	2.0	2-3	500	20
881 	Orange/Black	14 x 18	1.80	3-4	400	20


To Order: 800.511.6298 • www.medegenmed.com

Autoclavable Biohazardous Waste Bags

★ Indicates Made in the USA



844 (print)



5142.F (print)



844Y

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
5012	★	Red/Black	14 x 19	2.0	3-4	500	20
8-702	★	Orange/Black	19 x 23	2.0	7-10	200	20
882	★	Orange/Black	19 x 24	1.80	7-10	200	20
842	★	Red/Black	19 x 24	1.80	7-10	200	20
5013	★	Red/Black	19 x 24	2.0	7-10	200	20
8-150	★	Clear/Red	24 x 30	1.80	12-14	100	10
8-250	★	Clear/Red	24 x 36	1.80	13-16	100	10
854	★	Clear/Red	25 x 30	1.80	12-14	200	20
844	★	Red/Black	25 x 30	1.80	12-14	200	20
846	★	Red/Black	25 x 35	1.80	13-16	200	20
884	★	Orange/Black	25 x 35	1.80	13-16	200	20
8-200	★	Clear/Red	26 x 30	1.80	12-14	100	10
855	★	Clear/Red	31 x 38	1.80	20-30	100	10
856	★	Red/Black	31 x 38	1.80	20-30	100	10
858	★	Buff/Red	31 x 38	1.80	20-30	100	10
8-300	★	Clear/Red	36 x 35	1.80	33	100	10
5138	★	Buff/Red	38 x 46	2.0	44	100	10
5053	★	Clear/Red	38 x 47	1.75	44	100	10
872	★	Clear/Red	38 x 47	1.80	44	100	10
848	★	Red/Black	38 x 47	1.80	44	100	10
886	★	Orange/Black	38 x 47	1.80	44	100	10
864	★	Buff/Red	38 x 47	1.80	44	100	10
2120	★	Clear/Red	45 x 58	1.80	55	100	10
Max. Temp 285F Print/label: Biohazard Labeling, Biohazard Symbol, English/French.							
5020.F	★	Clear/Red	12 x 24	1.75	2-30	500	20
842Y	★	Yellow/Black	19 x 24	2.20	7-10	200	20
5142.F	★	Clear/Red	24 x 36	1.75	12-16	200	20
5025.F	★	Clear/Red	25 x 30	1.75	12-14	200	20
844Y	★	Yellow/Black	25 x 30	2.20	12-14	200	20
846Y*	★	Yellow/Black	25 x 35	2.20	13-16	200	20
856Y*	★	Yellow/Black	31 x 38	2.20	20-30	100	10
848Y*	★	Yellow/Black	38 x 47	2.20	44	100	10
No Print/Label Indicator,							
5023	★	Clear	38 x 47	1.80	44	100	10
5034	★	Clear	33 x 40	2.0	33	200	20

*Custom MTO, subject to MOQ.

Laundry and Linen Bags

A full line of Made in the USA printed and unprinted, color-coded Laundry and Linen Collection Bags provides a safe and efficient method of handling soiled or infectious linen and laundry from point of collection, through storage and transport, and into final processing.

- Wide range of standard sizes and film thickness, with and without print, provides maximum flexibility to meet your facility's needs. Custom Sizes and prints available. Contact us for a quote!
- Offered in both High Density and Linear Low Density resins for your specific application, bags are color-coded for easy identification by department
- Reinforced seams add durability and superior leak resistance
- Heavy Duty Linear Low Density Bags provide exceptional odor containment and are designed with superior tear and impact resistance for use in Multi-Story Laundry Chutes
- Typical areas of use: Housekeeping, Laundry, and most any treatment area or room where collection and transport of linens/clothing is required


★ Indicates **Made in the USA**

Laundry and Linen Bags - LLDPE Film, Star Seal, Coreless Roll


Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Soiled Linen.						
211M ★	White/Blue	25 x 34	1.20	12-16	250	25
C3102 ★	Blue/Black	30.5 x 41	.85	20-30	250	25
3102 ★	Blue/Black	30.5 x 41	1.00	20-30	250	25
285 ★	Blue/Blue	30.5 x 41	1.20	20-30	250	25
R285 ★	Light Blue/Blue	30.5 x 41	1.20	20-30	250	25
R285HA ★	Light Blue/Blue	30.5 x 41	1.05	20-30	250	25
264 ★	White/Blue	30.5 x 41	1.40	20-30	250	25
264M ★	White/Blue	30.5 x 41	1.20	20-30	250	25
3095 ★	Leaf Green/Black	40 x 46	1.00	40-45	200	20
51-40 ★	Blue/Black	40 x 46	1.20	40-45	100	10
Print/Label: Pediatric Soiled Linen.						
129M ★	Pink/Blue	30.5 x 41	1.10	20-30	250	25
Print/Label: Surgery Soiled Linen.						
203D ★	Mint Green/Black	30.5 x 41	1.40	20-30	250	25



264 (print)



129M (print)



203D (print)

Laundry and Linen Bags

★ Indicates Made in the USA



Laundry and Linen Bags - LLDPE Film, Star Seal, Coreless Roll

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Infectious Linen.						
208M ★	Yellow/Red	25 x 34	1.10	12-16	250	25
112M ★	Yellow/Red	30.5 x 41	1.10	20-30	250	25
112 ★	Yellow/Red	30.5 x 41	1.30	20-30	250	25
47-57 ★	Yellow/Red	31 x 41	1.00	20-30	250	25
47-07 ★	Yellow/Red	31 x 41	1.20	20-30	250	25
170M ★	Yellow/Red	38 x 45	1.10	44	250	25
51-45 ★	Yellow/Red	40 x 46	1.20	40-45	100	10



288

Laundry and Linen Bags - Draw Tape, LLDPE Film, Coreless Roll

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Soiled Linen.						
288 ★	Blue/Black	30.5 x 41	1.20	20-30	100	25



50-08

Laundry and Linen Bags - Draw Tape, LLDPE Film, Flat Pack

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
NO Print/Label.						
50-08	Blue	31 x 41	1.50	20-30	250	25



50-01

Anti-Static Kick Bucket Liners - Meet NFPA 56A, LLDPE Film, Star Seal, Coreless Roll

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Fits Oval and Round Kick Buckets. Print/Label: Anti-Static.						
50-01 ★	Blue/Black	18 x 8 x 16	1.30	5-7	500	50

Laundry and Linen Bags

★ Indicates Made in the USA


Laundry and Linen Bags - HDPE Film, Star Seal, Coreless Roll

Code#		Bag/Print Color	Size (inches)	Gauge/mic	Capacity/gal	Qty/cs	Qty/roll
Print/Label: Soiled Linen.							
RS304314PB	★	Blue/Black	30 x 43	14.00	20-30	250	25
RS304316PB	★	Blue/Black	30 x 43	16.00	20-30	250	25
3056	★	Blue/Black	40 x 46	14.00	40-45	250	25
Print/Label: Infectious Linen.							
RS304314Y	★	Yellow/Red	30.5 x 43	14.00	20-30	250	25
H355	★	Yellow/Red	30 x 54	14.00	30-35	200	20
RS384616Y	★	Yellow/Red	38 x 46	16.00	44	250	25



Unprinted Laundry and Linen Bags - LLDPE Film, Star Seal, Coreless Roll

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
NO Print/Label.							
301M	★	White/None	27 x 43	1.10	20-30	200	20
260M	★	Blue/None	30.5 x 41	1.20	20-30	250	25
R260M	★	Blue/None	30.5 x 41	1.20	20-30	250	25
201	★	Blue/None	30.5 x 41	1.40	20-30	250	25
264M	★	White/None	30.5 x 41	1.20	20-30	250	25
203	★	Mint Green/None	30.5 x 41	1.40	20-30	250	25
204XH	★	Clear/None	30.5 x 41	1.30	20-30	250	25
205M	★	Buff/None	30.5 x 41	.80	20-30	250	25
305	★	White/None	38 x 45	1.30	44	250	25
306	★	Blue/None	38 x 45	1.40	44	250	25
X1646	★	Mint Green/None	38 x 50	2.00	44	200	20



Unprinted Laundry and Linen Bags - HDPE Film, Coreless Roll

Code#		Bag/Print Color	Size (inches)	Gauge/mic	Capacity/gal	Qty/cs	Qty/roll
NO Print/Label.							
RS304314B	★	Blue/None	30.5 x 43	14.00	20-30	250	25
RS404816B	★	Blue/None	40 x 48	16.00	40-45	250	25

*Perforated bags on roll



RS304314B

Water Soluble Laundry Bags



Water Soluble Laundry Bags are designed for the safe and effective collection and handling of potentially contaminated or infectious linen and clothing. The bag is impervious to bacteria and viruses and is gas impermeable, exhibiting excellent odor containment properties. Contaminated linen is placed into the bag, sealed and placed directly into the washing machine, eliminating the need for staff to come in direct contact with the contents.

- Allows for easy isolation of contaminated and soiled linens to help reduce risk of cross-contamination
- 100% Biodegradable polyvinyl alcohol film dissolves completely into a non-toxic aqueous solution at 150°F during Hot Water machine washing and at energy saving 55°F for the Cold Water formulation
- Superior wet strength holdout capacity
- Integral tie strap for securing bag dissolves quickly to release contents during washing cycle
- Bags supplied on easy to dispense coreless rolls and individually packaged zip closure bags
- Sizes designed to fit most standard linen hampers
- Typical areas of use: in most any treatment area or room where potentially infectious linens or clothing come in contact with bodily fluids

Water Soluble Laundry Bags — Polyvinyl Alcohol Film, Coreless Roll

Code#	Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
Hot Water (150°F, 65°C), Heavy Duty. No Print/Label.						
1-344	Clear/None	19 x 22	1.00	7-10	100	25
1-341	Clear/None	26 x 33	1.00	20-25	100	25
1-347	Clear/None	28 x 39	1.00	30-32	100	25
1-348	Clear/None	36 x 39	1.00	40-45	100	25
Hot Water (150°F, 65°C), Standard Duty. No Print/Label.						
541-A	Clear/None	26 x 33	0.80	20-25	100	25
547-A	Clear/None	28 x 39	0.80	30-32	100	25
548-A	Clear/None	36 x 39	0.80	40-45	100	25
Cold Water (55°F, 13°C), Heavy Duty. No Print/Label.						
1-342	Clear/None	26 x 33	1.00	20-25	100	25
1-345	Clear/None	28 x 39	1.00	30-35	100	25
1-346	Clear/None	36 x 39	1.00	40-45	100	25

Hampers and Accessories

Hamper Stands facilitate convenient placement throughout the facility and ease of use. Stands are engineered for long-lasting performance in demanding health care environments with injection molded lid made from high impact ABS crack resistant resin to withstand repeated opening and closing.

- Hands free foot pedal actuation opens lid wide for easy access & closes automatically to help contain odors
- Solid 7/8" Chrome plated steel tubing with lower support platform for heavyweight applications
- Upper locking ring holds bags of at least 30 inches wide securely in place
- Open front design provides convenient access and bag removal
- Four 3" free turning casters — two with locking foot brakes — provide easy maneuverability
- Hamper pouch provides convenient storage for additional bags
- 19" wide x 20" deep x 37.5" tall. Available with or without labels



★ Indicates **Made in the USA**

Premium Chrome Hampers — Foot Pedal Activated Lid Operation

Code#	Description	Size (inches)	Capacity/gal	Qty/cs
15-9000A	Lid Label Pack: 5 label choices printed with Soiled Linen, Infectious Linen, Soiled Surgical Linen, Biohazardous Infectious Waste or Trash	19 x 20 x 37.5	30-33	1

Assorted Hardware and Accessories

Code#	Description	Capacity/gal	Qty/cs
15-24 ★	Toss-Away® Disposable Cardboard Hamper, yellow, no print, 17 x 25 x 15	20	10
15-9POUCH	Replacement Hamper Pouch, 14 1/2 x 14	1	
15-9LIDE	Premier Lid, Plastic, 19 x 21	1	

Hamper Labels (English, French, Spanish)

Code#	Print	Color
L411	Biohazardous Waste Infectious Waste	Red background w/White lettering
L412	Infectious Linen	Yellow background w/Red lettering
L413	Trash	Brown background w/White lettering
L414	Soiled Surgical Linen	Green background w/White lettering
L416	Soiled Linen	Blue background w/White lettering



Open Top Bags

Clear, LLDPE Film, Flat Pack

Minigrip Code #	Bag/Print Color	Size (inches) W x L or W x G x L	Qty/cs	Packaging
1.50mil Plain No Print				
XOPE10303	Clear/None	3 x 3	5,000	100
XOPE10304	Clear/None	3 x 4	5,000	100
XOPE10305	Clear/None	3 x 5	1,000	100
XOPE10406	Clear/None	4 x 6	1,000	100
XOPE10506	Clear/None	5 x 6	1,000	100
XOPE10508	Clear/None	5 x 8	1,000	100
XOPE15040208	Clear/None	4 x 2 x 8	1,000	100
XOPE10612	Clear/None	6 x 12	1,000	100
XOPE15040212	Clear/None	4 x 2 x 12	1,000	100
XOPE10709	Clear/None	7 x 9	1,000	100
XOPE10810	Clear/None	8 x 10	1,000	100
XOPE10912	Clear/None	9 x 12	1,000	100
XOPE11012	Clear/None	10 x 12	1,000	100
XOPE11216	Clear/None	12 x 16	1,000	100
XOPE11218	Clear/None	12 x 18	1,000	100
XOPE15080418	Clear/None	8 x 4 x 18	1,000	50
XOPE11420	Clear/None	14 x 20	250	25
XOPE11624	Clear/None	16 x 24	250	25
XOPE11824	Clear/None	18 x 24	250	25
XOPE12024	Clear/None	20 x 24	250	25
XOPE12430	Clear/None	24 x 30	250	25
2.0mil Plain No Print				
XOPE20304	Clear/None	3 x 4	5,000	100
XOPE20305	Clear/None	3 x 5	1,000	100
XOPE20406	Clear/None	4 x 6	1,000	100
XOPE20506	Clear/None	5 x 6	1,000	100
XOPE20508	Clear/None	5 x 8	1,000	100
XOPE25040208	Clear/None	4 x 2 x 8	1,000	100
XOPE20609	Clear/None	6 x 9	1,000	100
XOPE20612	Clear/None	6 x 12	1,000	100
XOPE20709	Clear/None	7 x 9	1,000	100
XOPE20810	Clear/None	8 x 10	1,000	100
XOPE20912	Clear/None	9 x 12	1,000	100
XOPE21012	Clear/None	10 x 12	1,000	100
XOPE21016	Clear/None	10 x 16	1,000	100
XOPE21216	Clear/None	12 x 16	1,000	50
XOPE21218	Clear/None	12 x 18	1,000	50
XOPE25080418	Clear/None	8 x 4 x 18	1,000	50
XOPE21420	Clear/None	14 x 20	250	25
XOPE21624	Clear/None	16 x 24	250	25
XOPE21824	Clear/None	18 x 24	250	25
XOPE22024	Clear/None	20 x 24	250	25
XOPE22430	Clear/None	24 x 30	250	25
4.0mil Plain No Print				
XOPE40303	Clear/None	3 x 3	1,000	100
XOPE40304	Clear/None	3 x 4	1,000	100
XOPE40305	Clear/None	3 x 5	1,000	100
XOPE40406	Clear/None	4 x 6	1,000	100
XOPE40508	Clear/None	5 x 8	1,000	100
XOPE40612	Clear/None	6 x 12	1,000	100
XOPE40709	Clear/None	7 x 9	1,000	100
XOPE40810	Clear/None	8 x 10	1,000	50
XOPE40912	Clear/None	9 x 12	500	50
XOPE41012	Clear/None	10 x 12	500	50
XOPE41016	Clear/None	10 x 16	500	50
XOPE41216	Clear/None	12 x 16	500	50
XOPE41218	Clear/None	12 x 18	500	25
XOPE45080418	Clear/None	8 x 4 x 18	500	25
XOPE41420	Clear/None	14 x 20	250	25
XOPE41624	Clear/None	16 x 24	250	25
XOPE41824	Clear/None	18 x 24	250	25
XOPE42024	Clear/None	20 x 24	250	10
XOPE42430	Clear/None	24 x 30	250	10

Simple easy to use open top bags. **No Closure** provides the flexibility to heat seal or use another closure style.



XOPE40912

Institutional Trash Can Liners

★ Indicates Made in the USA

Institutional Trash Can Liners provide a leak resistant containment solution for non-medical waste disposal and are used throughout the health care facility to collect and contain refuse. Available in Linear Low Density (LLDPE), ***TRANS4M™** PCR Content LLDPE and High density (HDPE) Polyethylene film.

- Available in a wide variety of standard sizes, film gauges, designed to meet your specific application. *Custom specifications available — contact us for a quote!*
- Made from premium resins with ultra-durable seams and puncture resistance.
- Clear liners allow quick and easy visualization of contents before disposal while colored options hide unsightly refuse
- LLDPE liners offer maximum strength and puncture/tear resistance due to their superior elasticity, provide excellent odor containment and are best suited for handling a variety of refuse including objects with straight corners/sharper edges
- **TRANS4M™** PCR Liners provide durable and trusted strength and containment performance AND are a blend of virgin LLDPE and post-consumer-recycled (PCR) content. Giving USED plastic bags and film a second life INSTEAD of sending them to landfills or the ocean.
- HDPE liners offer dependable durability and weight capacity in a lighter gauge film that is most cost effective for refuse without sharp or jagged edges
- Latex free and environmentally safe for incineration disposal
- Easy dispensing in convenient flat pack or coreless rolls
- Not for use with Biohazardous Waste



Institutional Trash Can Liners - **TRANS4M™** PCR + LLDPE Film, Star Sealed, Coreless Rolls

Transforming your plastic trash into plastic trash can liners!

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
PCR333909C	★	Clear/None	33 x 39	0.9	33	250	25
PCR333909K	★	Black / None	33 x 39	0.9	33	250	25
PCR384509C	★	Clear/None	38 x 45	0.9	44	250	25
PCR384509K	★	Black / None	38 x 45	0.9	44	250	25
PCR385809C	★	Clear/None	38 x 58	0.9	60	100	10
PCR385809K	★	Black / None	38 x 58	0.9	60	100	10
PCR434815C	★	Clear/None	43 x 48	1.5	56	100	10
PCR434815K	★	Black / None	43 x 48	1.5	56	100	10



Institutional Trash Can Liners

★ Indicates Made in the USA

Institutional Trash Can Liners – LLDPE Film, Star Sealed, Coreless Rolls

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs	Qty/roll
No Print or Label.							
23-01C	★	Clear	17 x 17	.35	3-4	1,000	50
4621	★	Brown	23 x 29	.50	8-10	500	50
4625	★	Clear	23 x 29	.50	8-10	500	50
4610	★	Brown	24 x 23	.45	8-10	1,000	50
4613	★	Clear	24 x 23	.45	8-10	1,000	50
4618	★	Brown	24 x 32	.40	12-16	1,000	50
4620	★	Clear	24 x 32	.40	12-16	1,000	50
T612	★	Clear	24 x 32	1.0	12-16	500	50
4627	★	Clear	30 x 36	.40	20-30	500	50
4637	★	Clear	30 x 36	.65	20-30	250	25
4630	★	Brown	30 x 36	.65	20-30	250	25
23-15W	★	White	30 x 36	.75	20-30	200	25
4644	★	Brown	33 x 39	.65	33	250	25
4646	★	Clear	33 x 39	.65	33	250	25
1290R	★	Clear	33 x 39	.90	33	250	25
T498	★	Clear	33 x 39	1.50	33	200	25
4648	★	Brown	33 x 39	.90	33	250	25
23-35C	★	Clear	33 x 40	1.25	33	100	10
T564	★	Brown	38 x 45	1.30	44	200	20
642	★	Clear	38 x 45	1.50	44	100	10
4661	★	Brown	38 x 55	.80	60	200	20
1212	★	Clear	38 x 58	.80	60	200	20
2750	★	White	38 x 58	1.25	60	100	20
R3906	★	Clear	40 x 46	.60	40-45	250	25
4662	★	Brown	40 x 46	.60	40-45	250	25
4660	★	Clear	40 x 46	.70	40-45	200	20
23-50C	★	Clear	40 x 46	1.10	40-45	100	10
3916	★	Black	40 x 48	1.50	40-45	200	20
3907	★	Brown	43 x 47	.80	55	100	10
4666	★	Clear	43 x 48	1.50	56	100	10
R4666	★	Clear	43 x 48	1.50	56	60	15
4665	★	Brown	43 x 48	1.50	56	100	10
Print: Recycle with Symbol							
X3005	★	Leaf Green	30.5 x 44	1.20	25-35	250	25
X3046	★	Leaf Green	43 x 47	1.20	55	250	25

Institutional Trash Can Liners – LLDPE Film, Flat Seal, Coreless Rolls (Custom-MTO)

Code#		Bag/Print Color	Size (inches)	Gauge/mil	Capacity/gal	Qty/cs
No Print/Label.						
T800	★	Clear	25 x 24 x 49	2.80	55-60	20

Institutional Trash Can Liners

★ Indicates Made in the USA

Institutional Trash Can Liners - HDPE Film, Star Seal, Coreless Interleaved Rolls

Code#	Bag/Print Color	Size (inches)	Gauge/mic	Capacity/gal	Qty/cs	Qty/roll
Print: Recycle with Symbol.						
RS242406N	Clear	24 x 24	6.00	8-10	1,000	50
RS242408N	Clear	24 x 24	8.00	8-10	1,000	50
RS243306N	Clear	24 x 33	6.00	12-16	1,000	50
RS243308N	Clear	24 x 33	8.00	12-16	1,000	50
RS303708N	Clear	30 x 37	8.00	20-30	500	25
RS303710N ★	Clear	30 x 37	10.00	20-30	500	25
RS303713N ★	Clear	30 x 37	13.00	20-30	500	25
RS303716N ★	Clear	30 x 37	16.00	20-30	500	25
RS334009N ★	Clear	33 x 40	9.00	33	500	25
RS334011N ★	Clear	33 x 40	11.00	33	500	25
RS334016N ★	Clear	33 x 40	16.00	33	250	25
RS334017N ★	Clear	33 x 40	17.00	33	250	25
RS334017W ★	White	33 x 40	17.00	33	250	25
RS366017N ★	Clear	36 x 60	17.00	55	200	25
RS374614N ★	Clear	37 x 46	14.00	44	250	25
RS374616N ★	Clear	37 x 46	16.00	44	250	25
RS386014N ★	Clear	38 x 60	14.00	60	200	25
RS386017N ★	Clear	38 x 60	17.00	60	200	25
RS404810N ★	Clear	40 x 48	10.00	40-45	250	25
RS404812N ★	Clear	40 x 48	12.00	40-45	250	25
RS404814N ★	Clear	40 x 48	14.00	40-45	250	25
RS404816N ★	Clear	40 x 48	16.00	40-45	250	25
RS404817N ★	Clear	40 x 48	17.00	40-45	250	25
RS404818W ★	White	40 x 48	18.00	40-45	200	25
RS404819K ★	Black	40 x 48	19.00	40-45	200	25
RS434814N ★	Clear	43 x 48	14.00	56	200	25
RS434816N ★	Clear	43 x 48	16.00	56	200	25

*Perforated bags on roll



Sharps Containers

Sharps-Tainer™

Our complete line of Stackable and Non-Stackable Sharps Containers is offered in a variety of sizes and configurations designed to give you the ability to choose the right product for your specific application. Regulated as Class II Medical devices and manufactured in accordance with FDA Quality System Regulations, as well as the International Quality System Standards ISO 13485:2003; all Medegen Medical Products Sharps Containers meet or exceed OSHA and Joint Commission Standards, CDC/NIOSH guidelines, ASTM requirements and DOT Regulations including 49CFR173.24, 49CFR173.24a, and 49CFR173.134.

- Puncture and Leak Resistant – Strong polypropylene construction meets or exceeds industry standards
- Locking Lids – Help keep contents secure for final disposal
- Autoclave/Incinerator Safe – Containers are Autoclave Safe up to a maximum of 280°F (137°C). When properly incinerated, containers emit only carbon dioxide and water
- Optional Wall Mounting Brackets and large container cart with locking casters
- Typical areas of use: Operating Room, Emergency Room, Same Day Surgery, Primary Care Treatment areas, Patient Room, Phlebotomy Labs, Chemistry & Hematology, Pathology, Toxicology, Microbiology, Immunization, Drug Screening, or wherever Biohazardous Sharps exist



Medegen's sharps containers are primarily made in the USA for your proven quality and dependable supply.

Look for the ★ next to items made in the USA.

Stackable Sharps Containers

★ Indicates **Made in the USA**

- Stackable design provides a convenient stable platform during use, economical storage, and lower shipping cost
- Wide horizontal drop lid models accommodate larger syringes and allow full utilization of fill line capacity
- Tortuous Path Lids help limit hand exposure for added security
- Transparent lids facilitate quick visual identification of fill level
- Translucent containers facilitate easy side viewing of content level
- Compact containers are conveniently sized for easy transport
- Extra large containers feature a wide hinged lid and round funneled vertical drop to accommodate endoscopic accessories, elongated disposable instruments and probes
- 5.4qt sharps fit into most 5 and 5.4qt sharps cabinets



Stackable Sharps-Tainer™ — Polypropylene

Code#		Container/Print Color	Size (inches)	Capacity	Qty/cs
Compact Print/Label: Biohazard, Biohazard Symbol.					
8701		Natural/Red/Black	6.5 x 2.5 dia.	.025 qt	48/cs
Small. Print/Label: Biohazard, Biohazard Symbol.					
8702	★	Red/Black	3 1/2 x 3 1/2 x 7	1 qt	72
8702T	★	Translucent Red/Black	3 1/2 x 3 1/2 x 7	1 qt	72
8702TY	★	Yellow/Black	3 1/2 x 3 1/2 x 7	1 qt	72
Medium. Print/Label: Biohazard, Biohazard Symbol.					
8703	★	Red/Black	10 x 7 x 5	1 gal	24
8703T	★	Translucent Red/Black	10 x 7 x 5	1 gal	24
Medium. Tortuous Path Lid. Print/Label: Biohazard, Biohazard Symbol.					
8706	★	Red/Black	10 x 7 x 6 3/4	1 gal	12
8706T	★	Translucent Red/Black	10 x 7 x 6 3/4	1 gal	12
Medium. Print/Label: Biohazard, Biohazard Symbol.					
8708	★	Red/Black	10 1/2 x 3 1/2 x 10 1/2	5.4qt	12
8708C	★	Clear/Black	10 1/2 x 3 1/2 x 10 1/2	5.4qt	12
8708T	★	Translucent Red/Black	10 1/2 x 3 1/2 x 10 1/2	5.4qt	12
8708TY	★	Yellow/Black	10 1/2 x 3 1/2 x 10 1/2	5.4qt	12
Patent Pending.					
Large. Print/Label: Biohazard, Biohazard Symbol.					
8704	★	Red/Black	10 x 7 x 9 1/2	2 gal	24
8704T	★	Translucent Red/Black	10 x 7 x 9 1/2	2 gal	24
8704TY	★	Yellow/Black	10 x 7 x 9 1/2	2 gal	24



Stackable Sharps Containers

★ Indicates Made in the USA

Stackable Sharps-Tainer™ — Polypropylene *(continued)*

Code#		Container/Print Color	Size (inches)	Capacity	Qty/cs
Large. Tortuous Path Lid. Print/Label: Biohazard, Biohazard Symbol.					
8707	★	Red/Black	10 x 7 x 11 ¹ / ₄	2 gal	10
8707T	★	Translucent Red/Black	10 x 7 x 11 ¹ / ₄	2 gal	10
8707TY	★	Yellow/Black	10 x 7 x 11 ¹ / ₄	2 gal	10
8710	★	Red/Black	12 ¹ / ₂ x 6 x 13 ¹ / ₂	3 gal	12
Large. Flat Lockable Sliding Lid. Print/Label: Biohazard, Biohazard Symbol.					
8710TF	★	Red/Black	12 ¹ / ₂ x 6 x 12	3 gal	12
Extra Large. Print/Label: Biohazard, Biohazard Symbol.					
8705	★	Red/Black	14 x 14 ³ / ₄ x 14 ¹ / ₄	8 gal	10
8705T	★	Translucent Red/Black	14 x 14 ³ / ₄ x 14 ¹ / ₄	8 gal	10
8716	★	Red/Black	14 x 14 ¹ / ₄ x 25 ³ / ₄	16 gal	6



Stackable Sharps-Tainer™ – Polypropylene, Chemo Waste Container

Chemo Waste (White) models also feature:

- Chemotherapy Waste Labeling/CYTA Symbol for easy identification
- Transparent top facilitates quick visualization of fill level

Code#		Container/Print Color	Size (inches)	Capacity/gal	Qty/cs
Extra Large. Print/Label: Biohazard, Biohazard Symbol.					
9751	★	White/Black & Yellow	14 x 14 ³ / ₄ x 14 ¹ / ₄	8	10
9752	★	White/Black & Yellow	14 x 14 ¹ / ₄ x 25 ³ / ₄	16	6



Medegen Stackable Sharps Family

Non-Stackable Sharps Containers

★ Indicates Made in the USA

Non-Stackable Sharps-Tainer™ – Polypropylene

Standard models feature:

- Funneled vertical drop and Extended Neck design help minimize needle exposure and overfilling
- Durable Polypropylene construction provides a stable base during use
- Extended handles help keep fingers away from opening during transport
- Needle removal port (#158, #178) provides convenient and compact disposal of sharps



Code#		Description	Color	Size (inches)	Capacity/qt	Qty/cs
Print/Label: Biohazard, Biohazard Symbol.						
158	★	Short Tray Size w/Needle Remover	Red/Black	3 1/2 x 3 1/2 x 5	0.7	40
185S	★	Short Tray Size, Locking Cap	Red/Black	3 1/2 x 3 1/2 x 5	0.7	40
178	★	Tall Tray Size w/Needle Remover	Red/Black	3 1/2 x 3 1/2 x 10	1.7	20
185	★	Tall Tray Size, Locking Cap	Red/Black	3 1/2 x 3 1/2 x 10	1.7	20
184	★	Small Size, Locking Cap	Red/Black	10 x 6 x 9	4.7	12
182	★	Medium Size, Locking Cap	Red/Black	10 x 6 x 11	6.2	12
180	★	Large Size, Locking Cap	Red/Black	10 x 6 x 15	10	12
184R	★	Rectangular, Locking Cap	Red/Black	10 1/4 x 6 x 7 3/4	4.7	12
186W	★	Large Wide Opening, Locking Cap	Red/Black	12 x 9 x 11	12	12
187	★	Extra Large Wide Opening, Locking Cap	Red/Black	12 x 9 x 18	23	4

Chemotherapy Waste Sharps-Tainer™

Chemo Waste (White) models feature:

- Chemotherapy Waste Labeling/ CYTA Symbol for easy identification

Code#		Description	Color	Size (inches)	Capacity/qt	Qty/cs
Print/Label: CYTA Symbol, Chemotherapy Waste.						
930	★	Small Size, Locking Cap	White/Black & Yellow	3 1/2 x 3 1/2 x 10	1.7	20
931W	★	Standard Size, Large Opening, Locking Cap	White/Black & Yellow	12 x 9 x 11	12	12
932	★	Large Size, Large Opening, Locking Cap	White/Black & Yellow	12 x 9 x 18	23	4



Non-Stackable Sharps Containers



Stackable Sharps Container Hardware — White

Code#	Description	Size (inches)	Qty/cs
8708BR	5.4QT. Sharps Wire Bracket	—	1
8780	Wall Mount for 8702 series	—	1
8782	Locking Wall Mount for 8703, 8704 series	—	1
8783	Locking Wall Mount for 8706, 8707 series	—	1
8784	Locking Wall Mount for 8705 series	—	1
8785	Locking Wall Mount for 8706, 8707 series	—	1
8786F	Locking Bracket for 3 Gallon 8710TF	—	1
8790	Cart w/Locking Casters for 8705 & 8716	13 ¹ / ₄ x 12 ¹ / ₂ x 15 ¹ / ₄	1

Non-Stackable Sharps Container Accessories

Code#	Description	Qty/cs
R175	Dovetail Bracket for use with 185, 178, 158, 185S	12



Solidifiers

★ Indicates Made in the USA

MEDEGEN MEDICAL PRODUCTS offers non-biodegradable solidifiers for various situations involving blood or other potentially infectious fluids. These solidifiers are formulated with a wetting agent which are ideal for encapsulating blood, infectious fluids found in suction canisters, kick buckets, and biohazardous bags.

- Solidifies liquids into an easy to handle gel
- Isolates, contains and prevent the spread of biohazardous fluids
- Minimizes splashing, spraying and splattering of blood or other potentially infectious fluids
- Eliminates the problem of disposing free liquids in landfills

RED Z

Code #		Description	Qty/cs
2030	★	Red-Z 8oz Bottle	6
2031	★	Red-Z Single Dose Bottle	75
2032	★	Red-Z Single Dose Bottle	75
P00-41105	★	Red-Z 4.4oz Needle Nose Bottle	24
P00-41102	★	Red-Z 11oz Shaker Bottle	12
2043M	★	Red-Z 16oz Bottle	15
P00-41175	★	Red-Z Tea Bag Pouch	100
P00-41131	★	Red-Z 1oz Diamond Pouch	100

GREEN Z

Code #		Description	Qty/cs
P00-42005	★	Green-Z Single Dose Bottle	288
P00-42010	★	Green-Z 5oz Shaker Bottle	24
M2042	★	Green-Z 22oz Bottle	12



Spill Control Solutions

★ Indicates **Made in the USA**

Emergency Response Kits comply with OSHA, CDC and State Health Department. Our Emergency Response Kits are a combination of Red-Z solidifiers and personal protection items for a safe, cost-effective clean-up of blood or body fluid spills. Encapsulating and containing a potentially infectious spills, helps reduce employee exposure and cross contamination to co-workers, patients, and their families. Find the kit that's right for you.

Bio-Hazard Spill Kits

Code #		Description	Qty/cs
2035	★	Kit: Biohaz Red-Z Emergency Response	36
2037	★	Kit: Biohaz Deluxe Emergency Response	6
2038	★	Kit: Biohaz Economy	72
2040	★	Kit: Bio Hazard Spill	16
P00-17100	★	Kit: Biohaz Universal Precaution	24
P00-17120	★	Kit: Biohaz EZ Cleans	24
P00-17121	★	Kit: Biohaz EZ Cleans Plus	24



Hazardous Drug Clean-up Kit

Code #		Description	Qty/cs
P00-49725	★	Hazardous Spill Clean-up	12



Spill Control Solutions

Biohazardous Spill Clean-Up Kits							Hazardous Drug Clean-Up Kit
	2035	2037	2038	P00-17100	P00-17121	P00-17120	P00-49725
Gown		Impervious		Apron			Chemo
Gloves	Vinyl	Vinyl	Vinyl	Vinyl	Vinyl	Vinyl	Nitrile
Face Mask		w/Safety Shield		w/Safety Shield			W/Safety Shield
Scoop & Scraper	Y	Y	Y	Y	Y	Y	Y
Absorbent Towel					Y	Y	
Liner(s)	Biohazardous Waste & LD Liner	Biohazardous Waste & LD Liner	Biohazardous Waste	Biohazardous Waste	Biohazardous Waste	Biohazardous Waste	Chemotherapy Waste
Absorbent Material	Red-Z	Red-Z	Red-Z	Red-Z	Red-Z	Red-Z	Green-Z
Wiper Pad	Y	Y					Y
Sanizide	Y	Y	Y	Y	Y		
P.A.W.S.	Y	Y	Y	Y	Y	Y	
Shoe Cover							Y

★ The above items are all **Made in the USA**

OSHA Compliance Guidelines

KITS	OSHA Compliance Guidelines
Item: 2037 Deluxe Emergency Response	29 CFR1910.1030 – Bloodborne Pathogens
Item: P00-17100 U.P. Compliance	

Protective Apparel



Gowns

Code #	Description	Qty/cs	Packaging
99916B	Isolation Gown, Blue, Poly	75	15



Vision-Tek™ Protective Eyewear

Code #	Description	Qty/cs	Packaging
*206-	Heavy Duty Safety Goggles, Clear	12	1
*208-	Safety Glasses, Clear	10	1
*209-	Safety Spectacles, Clear	36	1

*Meets the requirements of ANSI Z87.1

Product Number Index

109.....	25	2030.....	43	4914.....	7	117M.....	25	47-41.....	24
110.....	25	2031.....	43	4915.....	7	1290R.....	36	47-43.....	24
112.....	30	2032.....	43	4916.....	7	129M.....	29	47-45.....	24
116.....	25	2035.....	44	4921.....	24	1-341.....	32	47-55.....	25
117.....	25	2037.....	44	4923.....	24	1-342.....	32	47-57.....	30
119.....	25	2038.....	44	4925.....	7	1-344.....	32	47-71.....	23
158.....	41	2040.....	44	4930.....	24	1-345.....	32	47-73.....	25
169.....	25	2045.....	25	5001.4.....	27	1-346.....	32	47-CBL4.....	24
173.....	25	2120.....	28	5010.....	27	1-347.....	32	47-CBSI.....	18
178.....	41	2122.....	25	5012.....	28	1-348.....	32	48-07.....	20
180.....	41	2125.....	25	5013.....	28	15-24.....	33	48-10.....	20
182.....	41	2202.....	25	5023.....	28	15-9000A.....	33	4811B.....	9
184.....	41	2204.....	25	5034.....	28	15-9LIDE.....	33	4811G.....	9
185.....	41	2302.....	25	5051.....	19	15-9POUCH.....	33	48-19.....	20
187.....	41	2305.....	25	5053.....	28	170M.....	30	48-20.....	20
201.....	31	2634.....	23	5138.....	28	172M.....	25	48-29.....	20
203.....	31	2635.....	23	5139.....	21	184R.....	41	48-30.....	20
264.....	29	2750.....	36	6058.....	21	185S.....	41	48-39.....	20
285.....	29	3056.....	31	6059.....	21	186W.....	41	48-40.....	20
288.....	30	3095.....	29	6089.....	21	203D.....	29	4914.ORG.....	7
305.....	31	3102.....	29	6092.....	20	2043M.....	43	4915.ORG.....	7
306.....	31	3537.....	3	6220.....	20	204XH.....	31	49-90.....	8
642.....	36	3539.....	3	6456.....	21	205M.....	31	49-91.....	7
810.....	19	3549.....	3	6490.....	20	206.....	46	49-93.....	8
815.....	19	3660.....	6	6625.....	20	208.....	46	49-95.....	7
820.....	19	3670.....	7	8701.....	39	208M.....	30	49-96.....	8
825.....	19	3785.....	8	8702.....	39	209.....	46	49-99.....	7
830.....	19	3907.....	36	8703.....	39	211M.....	29	50-01.....	30
835.....	19	3916.....	36	8704.....	39	23-01C.....	36	50-02.....	26
840.....	19	4001.....	8	8705.....	40	23-15W.....	36	50-04.....	26
842.....	28	4002.....	8	8706.....	39	23-35C.....	36	50-08.....	30
844.....	28	4060.....	11	8707.....	40	23-50C.....	36	50-20.....	3
845.....	19	4075.....	7	8708.....	39	260M.....	31	5020.F.....	28
846.....	28	4210.....	24	8710.....	40	264M.....	29	50-21.....	4
848.....	28	4510.....	7	8716.....	40	264M.....	29	5025.F.....	28
852.....	27	4610.....	36	8780.....	42	301M.....	31	50-30.....	4
854.....	28	4613.....	36	8782.....	42	37-97.....	23	50-40.....	23
855.....	28	4618.....	36	8783.....	42	39-97.....	11	50-42.....	23
856.....	28	4620.....	36	8784.....	42	43-01.....	20	50-45.....	6
858.....	28	4621.....	36	8785.....	42	43-03.....	21	50-50.....	3
864.....	28	4625.....	36	8790.....	42	43-05.....	20	5050C.....	19
872.....	28	4627.....	36	9218.....	18	44-01.....	26	50-51.....	4
880.....	27	4630.....	36	9230.....	24	44-03.....	26	50-53.....	3
881.....	27	4637.....	36	9233.....	24	44-05.....	26	50-80.....	3
882.....	28	4644.....	36	9751.....	40	44-13.....	26	51-40.....	29
884.....	28	4646.....	36	9752.....	40	45-30.....	23	5142.F.....	28
886.....	28	4648.....	36	10-2011.....	26	45-43.....	23	51-45.....	30
930.....	41	4660.....	36	104M.....	25	45-45.....	23	52-04.....	19
932.....	41	4661.....	36	108M.....	25	45-50.....	23	541-A.....	32
1130.....	11	4662.....	36	112M.....	30	45-54.....	23	547-A.....	32
1212.....	36	4665.....	36	115H.....	25	45-75.....	24	548-A.....	32
2020.....	25	4666.....	36	117A.....	25	47-07.....	30	57-00.....	25

59-89	11	D4078	7	MGDZW270G	16	MGR4P1318	13	MGUV3P0406	17
59-93	11	D6624	20	MGDZW270J	16	MGR4P1320	13	MGUV3P0508	17
59-94	7	DSPB2020B	3	MGDZW270Q	16	MGR4P1424	13	MGUV3P0608	17
70-21	4	DSPB20C	3	MGR4P0202	12	MGR4PH02503	12	MGUV3P0808	17
70-31	4	DSPB20W	3	MGR4P0203	12	MGR4PH02506	12	MGUV3P0814	17
70-51	4	DTPB20W	4	MGR4P0205	12	MGR4PH0304	12	MGUV3P1212	17
8-150	28	F109	23	MGR4P0208	12	MGR4PH0404	12	MGZ2P0202	14
8-200	28	F110	23	MGR4P0303	12	MGR4PH0406	12	MGZ2P0203	14
8-250	28	F116	23	MGR4P0304	12	MGR4PH0408	12	MGZ2P0208	14
8-300	28	F117	23	MGR4P0305	12	MGR4PH0508	12	MGZ2P0303	14
842Y	28	F118	23	MGR4P0306	12	MGR4PH0604	12	MGZ2P0304	14
844Y	28	F124	23	MGR4P0404	12	MGR4PH0606	12	MGZ2P0305	14
846Y	28	F130	23	MGR4P0406	12	MGR4PH0608	12	MGZ2P0306	14
848Y	28	F135	23	MGR4P0407	12	MGR4PH1012	12	MGZ2P0404	14
856Y	28	F145	23	MGR4P0408	12	MGR4WH0304	13	MGZ2P0406	14
8-702	28	F149	23	MGR4P0505	12	MGR4WH0305	13	MGZ2P0407	14
8702T	39	F169	24	MGR4P0507	12	MGR4WH0406	13	MGZ2P0408	14
8702TY	39	F012THA*	5	MGR4P0508	12	MGR4WH0508	13	MGZ2P0505	14
8703T	39	F012WHI	5	MGR4P0510	12	MGR4WH0609	13	MGZ2P0507	14
8704T	39	F09WHI	5	MGR4P0604	12	MGR4WH0810	13	MGZ2P0508	14
8704TY	39	H355	31	MGR4P0606	12	MGR4WH0912	13	MGZ2P0510	14
8705T	40	HRD334014	26	MGR4P0608	12	MGR6P0406	13	MGZ2P0604	14
8706T	39	IP1215B3T	7	MGR4P0609	12	MGR6P0608	13	MGZ2P0606	14
8707T	40	IP2024B3T	7	MGR4P0610	12	MGR6P0810	13	MGZ2P0608	14
8707TY	40	IP46B3T	7	MGR4P0612	12	MGR6P0912	13	MGZ2P0609	14
8708BR	42	IP69B	7	MGR4P0808	12	MGR6P1012	13	MGZ2P0610	14
8708C	39	IP69BAS	8	MGR4P0810	12	MGR6P1215	13	MGZ2P0708	14
8708T	39	IP69BBL	9	MGR4P0812	12	MGR6P1318	13	MGZ2P0808	14
8708TY	39	IP69BG	9	MGR4P0912	12	MGS270G	16	MGZ2P0810	14
8710TF	40	IP69BY	9	MGR4P1012	12	MGS270J	16	MGZ2P0812	14
8786F	42	IP69RSTAT	8	MGR4P1013	12	MGS270L	16	MGZ2P0912	14
931W	41	IP69RSTAT2	8	MGR4P1212	12	MGS270Q	16	MGZ2P1010	14
95-CMS59	9	IP810RSTAT	8	MGR4P1215	12	MGST4P02503	17	MGZ2P1012	14
95-CMS810	9	L411	33	MGR4P1315	12	MGST4P0305	17	MGZ2P1013	14
99916B	46	L412	33	MGR4P1318	12	MGST4P0406	17	MGZ2P1212	14
AC1419C	27	L413	33	MGR4W0203	12	MGST4P0608	17	MGZ2P1215	14
AC1419R	27	L414	33	MGR4W0304	12	MGST4P0609	17	MGZ2P1315	14
AC1923C	27	L416	33	MGR4W0305	12	MGST4P0810	17	MGZ2P1318	15
AC1923R	27	M2042	43	MGR4W0406	12	MGST4P0912	17	MGZ2P1420	15
AC2530C	27	MGBD2P0202	13	MGR4W0508	12	MGST4P1212	17	MGZ2P1824	15
AC2530R	27	MGBD2P0203	13	MGR4W0609	12	MGST4P1318	17	MGZ2P2020	15
AC2535C	27	MGBD2P0304	13	MGR4W0810	12	MGSZW115S50	16	MGZ2W0203	15
AC2535R	27	MGBD2P0406	13	MGR4W0912	12	MGTZ1215	7	MGZ2W0304	15
AC3138C	27	MGBD2P0507	13	MGR4P0609	13	MGTZ69BLU	9	MGZ2W0305	15
AC3138R	27	MGBD2P0609	13	MGR4P0612	13	MGTZ69GRN	9	MGZ2W0406	15
AC3847C	27	MGBD2P0912	13	MGR4P0808	13	MGTZ69PUR	9	MGZ2W0508	15
AC3847R	27	MGCHEMO1215	18	MGR4P0810	13	MGTZ69RD	9	MGZ2W0609	15
BPL2214DT-W	21	MGCHEMO69	18	MGR4P0912	13	MGTZ69S	8	MGZ2W0810	15
C3102	29	MGCHEMO912	18	MGR4P1010	13	MGTZ69YEL	9	MGZ2W0912	15
CF2359	23	MGDZW175G	16	MGR4P1013	13	MGUV3P0203	17	MGZ2W1215	15
D2209	25	MGDZW175J	16	MGR4P1212	13	MGUV3P02509	17	MGZ4P02503	15
D2210	25	MGDZW175Q	16	MGR4P1215	13	MGUV3P0305	17	MGZ4P0304	15

Product Number Index

MGZ4P0305.....	15	RHPB20C.....	4	SBL2X912B.....	7	XOPE10912.....	34	Z2.0202.....	14
MGZ4P0404.....	15	RHPB20W.....	4	SG15WHI.....	5	XOPE11012.....	34	Z2.0203.....	14
MGZ4P0406.....	15	RS242406N.....	37	SG18WHI.....	5	XOPE11216.....	34	Z2.0203W.....	15
MGZ4P0508.....	15	RS242408N.....	37	SG9WHI.....	5	XOPE11218.....	34	Z2.0303.....	14
MGZ4P0606.....	15	RS242411RH.....	26	SPEC1-CC.....	10	XOPE11420.....	34	Z2.0304.....	14
MGZ4P0608.....	15	RS243306N.....	37	SPEC1-UV.....	11	XOPE11624.....	34	Z2.0304W.....	15
MGZ4P0609.....	15	RS243308N.....	37	SPEC-PNYB.....	10	XOPE11824.....	34	Z2.0305.....	14
MGZ4P0808.....	15	RS243311RH.....	26	SPEC-PNYG.....	10	XOPE12024.....	34	Z2.0305W.....	15
MGZ4P0810.....	15	RS303708N.....	37	SPEC-PNYP.....	10	XOPE12430.....	34	Z2.0404.....	14
MGZ4P0912.....	15	RS303710N.....	37	SPEC-PNYRT.....	10	XOPE15040208.....	34	Z2.0406.....	14
MGZ4P1010.....	15	RS303713N.....	37	SPEC-PNYY.....	10	XOPE15040212.....	34	Z2.0406W.....	15
MGZ4P1012.....	15	RS303716N.....	37	TI133.....	20	XOPE15080418.....	34	Z2.0408.....	14
MGZ4P1212.....	15	RS304314B.....	31	T498.....	36	XOPE20304.....	34	Z2.0507.....	14
MGZ4P1215.....	15	RS304314PB.....	31	T564.....	36	XOPE20305.....	34	Z2.0508.....	14
MGZ4P1318.....	15	RS304314RH.....	26	T612.....	36	XOPE20406.....	34	Z2.0508W.....	15
MGZ4P1424.....	15	RS304314Y.....	31	T800.....	36	XOPE20506.....	34	Z2.0606.....	14
MGZ4P1824.....	15	RS304316PB.....	31	TBPB20W.....	4	XOPE20508.....	34	Z2.0608.....	14
MGZ4P2020.....	15	RS334009N.....	37	THWIA.....	5	XOPE20609.....	34	Z2.0609.....	14
MGZ4W0305.....	15	RS334011N.....	37	UF95-1318BIOA.....	9	XOPE20612.....	34	Z2.0609W.....	15
MGZ4W0406.....	15	RS334016N.....	37	UF95-1517BIO.....	9	XOPE20709.....	34	Z2.0610.....	14
MGZ4W0508.....	15	RS334017N.....	37	UF95-600.....	11	XOPE20810.....	34	Z2.0808.....	14
MGZ4W0609.....	15	RS334017W.....	37	UF95-600BBL.....	10	XOPE20912.....	34	Z2.0810.....	14
MGZ4W0912.....	15	RS374614N.....	37	UF95-600BIO.....	9	XOPE21012.....	34	Z2.0912.....	14
MGZ6P1318.....	15	RS374616N.....	37	UF95-600BPM.....	9	XOPE21016.....	34	Z2.0912W.....	15
PCR3333909C.....	35	RS384616RH.....	26	UF95-600STB.....	10	XOPE21216.....	34	Z2.1012.....	14
PCR333909K.....	35	RS384816Y.....	31	UF95-600STBR.....	10	XOPE21218.....	34	Z2.1013.....	14
PCR384509C.....	35	RS386014N.....	37	UF95-610.....	11	XOPE21420.....	34	Z2.1212.....	14
PCR384509K.....	35	RS386017N.....	37	UF95-610BIO.....	9	XOPE21624.....	34	Z2.1215.....	14
PCR385809C.....	35	RS404810N.....	37	UF95-66SP.....	11	XOPE21824.....	34	Z2.1315.....	14
PCR385809K.....	35	RS404812N.....	37	UF95-66SPBIO.....	9	XOPE22024.....	34	Z2.1318.....	15
PCR434815C.....	35	RS404814N.....	37	UF95-BDES.....	10	XOPE22430.....	34	Z4.0304.....	15
PCR434815K.....	35	RS404816B.....	31	UF95-BL10.....	18	XOPE25040208.....	34	Z4.0406.....	15
P00-17100.....	44	RS404816N.....	37	UFDD1114.....	24	XOPE25080418.....	34	Z4.0508.....	15
P00-17120.....	44	RS404816RH.....	26	UFPE912.....	6	XOPE40303.....	34	Z4.0608.....	15
P00-17121.....	44	RS404817N.....	37	UFPSU1216.....	4	XOPE40304.....	34	Z4.0609.....	15
P00-41102.....	43	RS404818W.....	37	UFXC1012.....	19	XOPE40305.....	34	Z4.0808.....	15
P00-41105.....	43	RS404819K.....	37	UFXC1417.....	19	XOPE40406.....	34	Z4.0810.....	15
P00-41131.....	43	RS434814N.....	37	VGB-HPL.....	5	XOPE40508.....	34	Z4.0912.....	15
P00-41175.....	43	RS434816N.....	37	X1646.....	31	XOPE40612.....	34	Z4.1012.....	15
P00-42005.....	43	SBL2AP1215B.....	8	X2740.....	20	XOPE40709.....	34	Z4.1212.....	15
P00-42010.....	43	SBL2AP69B.....	8	X3005.....	36	XOPE40810.....	34	Z4.1215.....	15
P00-49725.....	44	SBL2AP810B.....	8	X3009B.....	24	XOPE40912.....	34	Z4.1318.....	15
R175.....	42	SBL2X1118B.....	7	X3046.....	36	XOPE41012.....	34	ZLAB120.....	7
R260M.....	31	SBL2X1215B.....	7	XOPE10303.....	34	XOPE41016.....	34	ZLABTRANSPORT.....	7
R285.....	29	SBL2X46B.....	7	XOPE10304.....	34	XOPE41216.....	34		
R285HA.....	29	SBL2X69.....	8	XOPE10305.....	34	XOPE41218.....	34		
R3906.....	36	SBL2X69B.....	7	XOPE10406.....	34	XOPE41420.....	34		
R3S66017N.....	37	SBL2X69B5.....	7	XOPE10506.....	34	XOPE41624.....	34		
R4666.....	36	SBL2X69S.....	8	XOPE10508.....	34	XOPE41824.....	34		
RD610.....	24	SBL2X810B.....	7	XOPE10612.....	34	XOPE42024.....	34		
RD630.....	24	SBL2X810B5.....	7	XOPE10709.....	34	XOPE42430.....	34		
RD650.....	24	SBL2X88B.....	7	XOPE10810.....	34	XOPE45080418.....	34		



MEASUREMENT & COLLECTION

PATIENT BEDSIDE ESSENTIALS

COLLECTION & CONTAINMENT

209 Medegen Drive
Gallaway, TN 38036
800.511.6298

Customer Service Fax:
631.714.6311

Customer Service:
CSservice@medegenmed.com

www.medegenmed.com



INTEPLAST GROUP
Healthcare



*To view our products online, please visit us at
www.medegenmed.com or directly with your
device by scanning the QR code to the left.*