



Specimen  
**TRANSPORT BAGS**

# Specimen Transport Bags

Minigrip / Medegen offers a comprehensive range of high-quality specimen transport bags designed to support the safe collection, storage, identification, and transport of laboratory specimens and cultures. These bags are ideal for use in hospitals, surgery centers, laboratories, and off-site facilities.

Each bag features a dual-pouch system that effectively separates and protects specimens from accompanying documentation, ensuring safety, organization, and compliance during handling and transport.

## Trusted Brands

Minigrip specimen packaging includes the trusted brands:

**Lab Guard • Speci-Zip • SAFEGUARD**

These brands have been relied upon by hospitals and laboratories for over 25 years.

## Features & Benefits

- FDA registered for quality and compliance
- Tear-Zone EZ-Open feature for quick and easy access
- Expanded temperature and STAT selector checkboxes for enhanced versatility
- Available in multiple standard sizes
- Options include clear or colored film
- Custom size, film color, print available for your unique application

## Closure types:

- Easy-to-use zip closure
- Reclosable adhesive
- Permanent tamper-evident adhesive

## Applications

Designed for a wide range of clinical and laboratory environments, including:

- Phlebotomy labs
- Chemistry and hematology
- Pathology
- Toxicology
- Microbiology
- Emergency rooms
- Treatment areas
- Drug screening
- Routine clinical specimen collection



4510



4914



SBL2X69B5

## Lab Guard® Speci-Zip® & SAFEGUARD® Lab Specimen Transport Bags — Zip Closure

Lab Guard®, Speci-Zip® & SAFEGUARD® brands provide convenient, easy-to-use, cost-saving features you just can't get with any other bag in a range of combinations. Choose below from single zip and liquid tight double zip closures, standard or cost saving destroyable biohazard symbol, TearZone® EZ open feature or NO TearZone, paperwork pouch flap or NO flap, STAT+Temp or just Temp selector boxes and more!

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SBL2X46B		4x6	Single Zip	Clear/Blk & Org		X	Standard		no flap	DispBag	20/50	1000
IP46B3T		4x6	Single Zip	Clear/Blk & Org		X	Standard	Score	w/flap	Zip Bag	10/100	1000
	3670	4x6	Single Zip	Clear/Blk & Org	X		Standard		w/flap	Zip Bag	10/100	1000
	49-91	6x6	Single Zip	Clear/Blk & Red	X		Standard		w/flap	Zip Bag	10/100	1000
SBL2X69B5		6x9	Double Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
SBL2X69B		6x9	Double Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBOX	2/500	1000
IP69B		6x9	Single Zip	Clear/Blk & Red			Standard	Score	w/flap	Zip Bag	10/100	1000
	4075	6x9	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4914	6x9	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4914.ORG	6x9	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	no flap	Zip Bag	10/100	1000
	D4078	6x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4510	7x9	Slider Zip	Clear/Blk & Org	X		Standard		w/flap	Zip Bag	10/50	500
SBL2X88B		8x8	Single Zip	Clear/Blk & Org		X	Standard		no flap	DispBag	20/50	1000
	49-95	8x8	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X810B5		8x10	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
SBL2X810B		8x10	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBOX	2/500	1000
	4915	8x10	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4915.ORG	8x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X912B		9x12	Single Zip	Clear/Blk & Org		X	Standard	TearZone®	no flap	DispBag	20/50	1000
	49-99	9x12	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4925	10x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X1118B		11X18	Single Zip	Clear/Blk & Org		X	Standard		no flap	Zip Bag	20/50	1000
SBL2X1215B		12x15	Single Zip	Clear/Blk & Org		X	Standard		no flap	Zip Bag	40/25	1000
MGTZ1215		12x15	Single Zip	Clear/Blk & Org		X	Standard	TearZone®	no flap	DispBag	5/100	500
IP1215B3T		12x15	Single Zip	Clear/Blk & Org		X	Standard	Score	w/flap	Zip Bag		250
	4916	12x15	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	5/100	500
	59-94	12x13x4	Single Zip	Clear/Blk & Red	NONE	NONE	Standard		w/flap	Zip Bag	10/100	1000
ZLABTRANSPT		12.75x12x6	Single Zip	Clear/Blk & Org	NONE	NONE	Standard		no doc pch	Zip Bag	5/100	500
ZLAB1420		14x20	Single Zip	Clear/Blk & Org	NONE	NONE	Standard		no doc pch	Zip Bag	5/100	500
IP2024B3T		20x24	Single Zip	Clear/Blk & Org		X	Standard		w/flap	Zip Bag	1/100	100

## Lab Guard® & Speci-Zip® Lab Specimen Transport Bags — Zip Closure, STAT

Lab Guard® & Speci-Zip STAT Specimen Bags are available in yellow, red, and blue color schemes. Other available features include cost-saving TearZone®, destroyable biohazard symbol or standard biohazard symbol and DispenserBag® or standard zip bag packaging.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
	3785	6x9	Single Zip	Yel Opaque/Blk		X	Standard	TearZone®	w/flap	Zip Bag	10/100	1000
MGTZ69S		6x9	Single Zip	Red Tint/Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
IP69RSTAT2		6x9	Single Zip	Red Tint/Blk		X	Standard	Score	w/flap	Zip Bag	10/100	1000
	4001	6x9	Single Zip	Red Tint/Blk		X	Standard	TearZone®	w/flap	Zip Bag	10/100	1000
IP69RSTAT		6x9	Single Zip	Red Opaque/Blk		X	Standard	Score	w/flap	Zip Bag	10/100	1000
SBL2X69S		6x9	Single Zip	Clr/Red,Wht,Blk		X	Standard		no flap	DispBag	20/50	1000
	4002	6x9	Single Zip	Blu Opaque/Blk		X	Standard	TearZone®	w/flap	Zip Bag	10/100	1000
IP810RSTAT		8x10	Single Zip	Red Opaque/Blk		X	Standard	Score	w/flap	Zip Bag	10/100	1000



MGTZ69S



SBL2X69S



4001



3785



4002



SBL2AP69B

## Lab Guard® & Speci-Zip® Lab Specimen Transport Bags — Zip Closure with Absorbent Pad Insert (3x5 Pad)

These Lab Guard® & Speci-Zip Specimen Bags feature a super-absorbent 3x5 absorbent pad to soak up leaking specimen cup fluids.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SBL2AP69B		6x9	Double Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
IP69BAS		6x9	Single Zip	Clear/Blk & Red			Standard	Score	w/flap	Zip Bag	10/100	1000
SBL2AP810B		8x10	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
SBL2AP1215B		12x15	Single Zip	Clear/Blk & Org		X	Standard		no flap	DispBag	20/50	1000

49-90



SBL2X69



## Lab Guard® Lab Specimen Transport Bags — Zip Closure, Utility (No Biohazard Print)

Lab Guard® UTILITY Specimen Bags are not printed with a biohazard symbol for hospitals and labs with in-facility processing and no transportation outside of the facility.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
	49-90	6x6	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000
SBL2X69		6x9	Double Zip	Clear/Blk & Org	NONE	NONE	NONE	TearZone®	no flap	DispBag	20/50	1000
	49-96	6x9	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000
	49-93	8x8	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000

# Lab Guard® & Speci-Zip® Lab Specimen Transport Bags — Zip Closure, Color-Codes

Lab Guard® & Speci-Zip® Color-Code Specimen Bags feature a range of color blocks on the front of bag to help hospitals and labs segregate specimens with specific processing needs. Features cost-saving TearZone®, destroyable biohazard symbol.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
MGTZ69BLU	4811B	6x9	Single Zip	Clear/Blu,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
IP69BBL		6x9	Single Zip	Blue/Red,Blk			Standard	Score	w/flap	Zip Bag	10/100	1000
MGTZ69YEL		6x9	Single Zip	Clear/Yel,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
IP69BY		6x9	Single Zip	Yellow/Red,Blk			Standard	Score	w/flap	Zip Bag	10/100	1000
MGTZ69GRN	4811G	6x9	Single Zip	Clear/Grn,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
IP69BG		6x9	Single Zip	Green/Red,Blk			Standard	Score	w/flap	Zip Bag	10/100	1000
MGTZ69RD		6x9	Single Zip	Clear/Red,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
MGTZ69PUR		6x9	Single Zip	Clear/Pur,Org,Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000



MGTZ69RED



MGTZ69BLU




MGTZ69YEL



MGTZ69GRN



MGTZ69PUR



## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure

Speci-Gard® Specimen Bags feature the innovative, easy-to-use, liquid-tight PRESS AND CLOSE® permanent adhesive closure, TearZone® EZ open system, DispenserBag® convenience packaging, and in some cases, a cost-saving destroyable biohazard symbol.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
95-CMS59		5x9	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-66SPBIO		6x6	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	no pouch	DispBag	100	2000
UF95-600BIO		6x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-600BPM		6x10	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	DispBag	100	1000
95-CMS810		8x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-610BIO		10x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-1318BIOA		13x18	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	Zip Bag	50	500
UF95-1517BIO		15x17	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	Zip Bag	50	500



UF95-600BIO



UF95-600BPM



UF95-600STB



UF95-600STBR

To Order: 800.533.1931  
www.medegenmed.com



## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure - STAT

Speci-Gard® STAT Specimen Bags are available in red and yellow film colors with all featuring the easy-to-use, liquid-tight, PRESS AND CLOSE® permanent adhesive closure, TearZone® EZ open system, cost saving destroyable biohazard symbol and DispenserBag® convenience packaging.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SPEC-PNYRT		6x10	Adhesive	Red Tint/Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
UF95-600STBR		6x10	Adhesive	Red Opaque/Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
UF95-600STB		6x10	Adhesive	Yel Opaque/Blk		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure & Absorbent Pad Insert (3x5 Pad)

Speci-Gard® Specimen Bags with absorbent pads are available with either single adhesive pouch (UF95-BDES) or DUAL-adhesive pouch (SPECI-CC). Item SPECI-CC secures both the specimen and the paperwork in separate pouches so that neither can be tampered with, without destroying the bag (chain-of-custody).

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
UF95-BDES		6x10	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	DispBag	10/100	1000
SPECI-CC		6x10	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	Adhesive	DispBag	10/100	1000



UF95-BDES



SPEC-PNYB



SPEC-PNYG



SPEC-PNYP

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure - Color Codes

Speci-Gard® Color-Code PRESS AND CLOSE® Adhesive Closure Specimen Bags feature a range of film colors to help hospitals and labs segregate specimens with specific processing needs.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SPEC-PNYB		6x10	Adhesive	Blu Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
UF95-600BBL		6x10	Adhesive	Blu Opaq/Blk & Red		X	Standard	TearZone®	w/flap	DispBag	10/100	1000
SPEC-PNYG		6x10	Adhesive	Grn Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
SPEC-PNYP		6x10	Adhesive	Pur Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000
SPEC-PNYY		6x10	Adhesive	Yel Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	10/100	1000



SPEC-PNYY



UF95-600BBL



UF95-600

To view our products online, please visit us at [www.medegenmed.com](http://www.medegenmed.com) or directly with your device by scanning the QR code below.



## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure, Utility (No Biohazard Print)

Speci-Gard® UTILITY Specimen Bags are not printed with a biohazard symbol for hospitals and labs with in-facility processing and no transportation outside of the facility. All bags offer TearZone® EZ-Open Feature and DispenserBag® Convenience Packaging.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
UF95-66SP		6x6	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	10/100	2000
UF95-600		6x10	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	10/100	1000
UF95-610		10x10	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	10/100	1000

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure, UV Protection

Speci-Gard® style PRESS AND CLOSE® Permanent Adhesive Closure Bag features opaque gray white film that prevents light from entering the bag, protecting light sensitive contents. Access to contents is available with EZ-Tear Perforation.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SPECI-UV		6x10	Adhesive	Wht Opq/Blk & Red	NONE	NONE	Standard	Perforation	w/flap	DispBag	10/100	1000

## Lab Specimen Transport Bags: Permanent Adhesive Closure, Foldover Closure

Lab specimen bags feature a traditional foldover permanent adhesive closure and score line for opening during processing.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
	59-89	6x9	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	10/100	1000
	4060	6x10	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	10/100	1000
	59-93	10x10	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	10/100	1000

## Lab Specimen Transport Bags: Adhesive Closure, Specialty Applications

Item 1130 features rugged welded side seams and permanent adhesive closure for pneumatic tube systems.

Item 39-97 features a reclosable adhesive closure for worry-free post-closure access.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
	1130	6x9	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Tear Notch	no flap	Bag	10/100	1000
	39-97	6x10	Reclosable Adh	Clear/Blk & Red	NONE	NONE	Standard	Pull Open	no flap	Zip Bag	10/100	1000



SPECI-UV



4060



1130



39-97



360 Motor Pkwy, Suite 800 • Hauppauge, NY 11788 • 800.533.1931  
 Customer Service Fax: 631.714.6311 • Customer Service: mgservice@medegenmed.com  
 ©2026 Medegen Medical Products

[www.medegenmed.com](http://www.medegenmed.com)

LIT# LIT0122

