



HOSPITAL FURNITURE

*Quality, Comfort and
Convenience*



HOSPITAL FURNITURE

PATIENT BEDSIDE CABINETS

Our patient Bedside Cabinets come assembled and are offered in ABS plastic.

- Color offered in ABS: Blue and Brown
- ABS plastic Bedside Cabinet top has a lip on 3 sides preventing item from falling off
- Left or Right mountable hanger for either Bedside Cabinet
- ABS plastic comes with a slide out shelf, drawer, closet and a vented shelf
- 4 x 2" casters, front casters are locking
- Dimension: L19" x W 19" x H 34"



PATIENT OVERBED TABLE

Medegen's Overbed Tables are designed to be:

- STABLE – with an “H” shaped base and four casters (2 x 1" and 2 x 2")
- DURABLE – we offer two options: wood grain laminate or ABS Plastic
- FUNCTIONAL – our Overbed Table has a gas spring lifting column with a fingertip self-locking handle enabling Healthcare workers to adjust the table height to accommodate the Patient Bed, Stretcher, Chair/Recliner or a Wheelchair.
ABS overbed table adjustable height: 31"~43"
Wood Laminate overbed table adjustable height: 30"~41"
- Overbed Table Dimensions:
ABS: L33" x W18" x H31"~43"
Wood Laminate: L35" x W16" x H30"~41"
- Overbed Tables are available in ABS plastic color Blue and Wood Laminate color Light Maple
- Table Weight capacity is 60 LBS



RECLINER-SLEEPER

5 Operating Handles

Operating Handles are on the left and right side of the recliner, making it easy to move or assisting with converting the recliner to a sleeper.



2 Safety Knobs

Safety Knobs are located on the left and right sides of the recliner-sleeper. Either knob can easily convert the recliner to a sleeper and back to a recliner.




3 Sleeper

Recliner-Sleeper has an ergonomic sleeping surface with drop down armrests providing an unobstructed sleeping area.




4 Gas Spring Back

The recliner has a Gas Spring back making the transition from recliner to sleeping surface effortless.



1 Arm Rests

Arm rests are made of soft comfortable foam.



6 Braking Caster

Recliner comes with 4 x 3" casters, 2 are locking and are toe operated.



Medegen Recliner-Sleepers are offered in 3 different sizes:

SMALL*

Chair: L 34" x W 24" x H 39"

Bed: L 77" x W 24" x H 18"

Sleeping surface: L 73" x W 21"

MEDIUM*

Chair: L 34" x W 27" x H 39"

Bed: L 79" x W 27" x H 18"

Sleeping surface: L 75" x W 22"

LARGE

Chair: L 35" x W 30" x H 42"

Bed: L 85" x W 30" x H 18"

Sleeping surface: L 80" x W 24"

BEIGE

*Available by special order

To Order: 800.511.6298

www.medegenmed.com



FDA



OTHER GREAT PRODUCTS FROM MEDEGEN


the **KALIKART**™ series of **HEALTHCARE CARTS**



See our full line of carts in this
informative trifold product brochure!



CHEVALIER **Hydraulic STRETCHER**



See all the features in our
comprehensive trifold
product brochure!



To view our products online, please visit us at www.medegenmed.com
or directly with your
device by scanning the QR code to the right.



360 Motor Pkwy, Suite 800 • Hauppauge, NY 11788 • 800.511.6298
Customer Service Fax: 631.714.6311 • Customer Service: CSservice@medegenmed.com
International Customer Service: IntCSservice@medegenmed.com

www.medegenmed.com

