

Biological COLLECTION

Starswab II™ collection is a complete bacteriology and transport system. The color coded cap has a protective over skirt which helps ensure the containment of biohazardous material.

- All Starswab II™ packaging is labeled with product number, lot # and expiration date.
- The vials are made of polypropylene with a thick tube wall to protect medium from dehydration and the spiral media chamber keeps sample immersed and protected during transport.
- The round bottom tube makes for easy racking and the seamless tube molding ensures tube strength and integrity.
- The rayon tip is inert, non-toxic and permits good sample retrieval and absorption.
- Amies Clear provides excellent survival for a wide range of bacteria and is compatible with gram stain.
- Amies Clear provides excellent environment for the survival of anaerobes through the complete immersion of the swabs into the gel and by the presence of the reducing agent, sodium thioglycollate
- Modified Liquid Stuart's and Liquid Amies provide excellent recovery in routine bacteriology and are compatible with gram stain and many rapid screening methods.

Starswab II™ packaging promotes optimal medium performance. A premium thick bag creates a barrier environment to prevent premature medium dehydration. The metalized bag further protects the medium from harmful UV rays. The bag also features a ziplock closure which simplifies resealing to better maintain the remaining swabs at peak performance.

Code #	Medium	Swab	Cap Color	Sterile	Case Pack	Packaging
P02-SP130X	Amies Clear	Single, Plastics Shaft, Rayon	Magenta	Y	500	10 Bags of 50 each
P02-SP130D	Amies Clear	Double, Plastic Shaft, Rayon	White	Y	500	10 Bags of 50 each
*P02-SP130L-D	Amies Clear w/sponge	Double, Plastic Shaft, Rayon	White	Y	500	10 Bags of 50 each
P02-SP132	Liquid Stuart's	Single, Plastic Shaft, Rayon	Teal	Y	500	10 Bags of 50 each
*P02-SP132D	Liquid Stuart's	Double, Plastic Shaft, Rayon	Yellow	Y	500	10 Bags of 50 each
*P02-SP132-W	Liquid Stuart's	Single, Straight Wire Shaft, Fine Tip	Orange	Y	500	10 Bags of 50 each

*Item is new Medegen Medical Product offering

MEDEGEN
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036
Toll Free 800-511-6298
Customer Service Fax 631-714-6311
E-mail cservice@medegenmed.com



INTEPLAST GROUP
Inteplast Healthcare



P02-SP130L-D



P02-SP132-W

